5266041205

Self-Adhesive Composite Ionomer Luting Cement

TOKUYAMA IONOTITE F

REF RC04

ENGLISH

Read all information, precautions and notes before using.

■ **PRODUCT DESCRIPTION AND GENERAL INFORMATION**

- TOKUYAMA IONOTITE F is self cured, radiopaque, fluoride releasing, self-adhesive composite ionomer luting cement system that has excellent adhesive properties toward precious metal alloys, non-precious metal alloys and teeth without requiring separate etching, priming or bonding steps.
- The main components of the TOKUYAMA IONOTITE F Kit are powder (the "Powder") and liquid (the "Liquid") that form a white cement (the "Cement") when combined. The Powder contains Fluoroaluminosilicate glass filler, Silica filler and Benzyl peroxide (BPO). The Liquid contains 2-Hydroxyethyl methacrylate (HEMA), 1,6-bis(methacryloxyhexyloxy)trimethyl hexane (UDMA), Phosphoric acid monomer and 6-Methacryloyloxyhexyl 2-thiouracil-5-carboxylate (MTU-6).

■ **INDICATIONS**

- TOKUYAMA IONOTITE F is indicated for permanent cementation of:
 - porcelain fused to metal (PFM) crowns and bridges,
 - cast crowns and bridges with resin facings,
 - metal inlays, onlays, crowns and bridges,
 - cast posts.
- TOKUYAMA IONOTITE F is **NOT indicated** for:
 - cementation of all ceramic crowns, all resin inlays or porcelain veneers because TOKUYAMA IONOTITE F has linear expansion of approximately 2% that may lead to fracture of these restorations,
 - cementation of prefabricated post, core build-up, restorative fillings or fixation of loose teeth because oxygen located at the cement surface inhibits hardening of TOKUYAMA IONOTITE F,
 - bonding of Maryland bridges or orthodontic bracket/retainer.

■ **CONTRAINDICATIONS**

TOKUYAMA IONOTITE F contains methacrylic monomers. TOKUYAMA IONOTITE F is contraindicated for patients allergic to or hypersensitive to methacrylic and related monomers and should NOT be used for such patients.

■ **PRECAUTIONS**

- DO NOT use TOKUYAMA IONOTITE F for any purpose other than those listed in these instructions. Use TOKUYAMA IONOTITE F only as directed herein.
- TOKUYAMA IONOTITE F is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
- DO NOT use TOKUYAMA IONOTITE F if the safety seals are broken or appear to have been tampered with.
- IF TOKUYAMA IONOTITE F causes an allergic reaction or oversensitivity, discontinue its use immediately.
- The use at all times of examination gloves (plastic, vinyl or latex) is recommended to avoid the possibility of allergic reactions from methacrylic monomers. Certain ingredients may penetrate through examination gloves. If a portion of TOKUYAMA IONOTITE F comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- Avoid contact of TOKUYAMA IONOTITE F with eyes, skin and clothing. Minimize contact with mucosal membrane.
 - IF TOKUYAMA IONOTITE F comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - IF TOKUYAMA IONOTITE F comes into contact with the skin or clothing, immediately saturate the area with an alcohol-soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
- TOKUYAMA IONOTITE F should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- To avoid being swallowed by mistake, do not leave TOKUYAMA IONOTITE F unsupervised within the reach of patients and children.
- Clean the plastic spatula and spoon with alcohol after use.
- Do not expose Liquid or its vapor to open flame.

■ **INDICATIONS FOR WORKING TIME & SETTING TIME**

- Mixing of Powder with Liquid should be completed within 30 seconds to maximize the working time.
- Working time depends on the Powder/Liquid ratio and temperature. The following table illustrates the relationship between the Powder/Liquid ratio and working time at 23°C(73°F). Note: Warm temperatures shorten the working time.

Relationship between the Powder/Liquid ratio and working time:

The Powder/Liquid Ratio (g/g)	2.6	2.8 (standard ratio)	3.0
Working time at 23°C (73°F)	3min00sec	2min20sec	2min00sec

- To keep the film thickness thin, seat the prosthesis as soon as possible after the application of cement to the tooth surface because tooth temperature accelerates the setting time. Special attention is required when seating a long-span bridge, a cast post or multiple restorations at one sitting.
- Waxy excess cement can be easily removed within 1 to 3 minutes after seating. The ideal time is 2 minutes. After 4 minutes the cement will take considerable effort to remove.
- Maintain normal occlusal pressure on the prosthesis with a cotton roll for 6 or more minutes from initial placement to complete set.
- Sufficient adhesion to enamel occurs 10 or more minutes from initial placement. Adhesion to dentin occurs 20 or more minutes from initial placement.

■ **PRECAUTIONS FOR MEDICATIONS AND MATERIALS**

- Some materials and medicaments (hemostatic agent) inhibit setting/adhesion of TOKUYAMA IONOTITE F for an extended period even after careful meticulous cleansing with water. DO NOT USE products which contain:
 - eugenol,
 - iodine (iodine tincture),
 - diammine silver fluoride [molecular formula: Ag(NH₃)₂F],
 - phenols such as parachlorophenol, guaiacol, phenol,
 - aluminum chloride,
 - ferric sulfate.
- TOKUYAMA IONOTITE F does not adhere to tooth structure immediately after applications of materials which contain:
 - hydrogen peroxide (oxydol),
 - sodium hypochlorite,
 - chlorhexidine gluconate.

We advise avoiding the use of these materials on the same day of cementation. However these materials can be used for root canal treatment because the inhibition caused by these materials usually disappear within five days.

- TOKUYAMA IONOTITE F does not adhere to tooth structure in case the bonding surface is contaminated with medicaments (hemostatic agent) which contain:
 - aluminum sulfate,
 - epinephrine.

DO NOT apply these medicaments to the bonding surface. Thoroughly rinse the bonding surface with water if the bonding surface is contaminated.

■ **PROCEDURAL GUIDELINES FOR CAST POST**

- Root canal must be dried with absorbent points before cementing because excess water in the canal may cause premature detachment.
- DO NOT prepare a cast post abutment immediately after setting because insufficient adhesion may cause premature detachment. In order to achieve sufficient adhesion to dentin, maintain appropriate pressure onto the cast post for 6 minutes and keep for 14 minutes without pressure (sum total 20 minutes is required).

■ **STORAGE**

- Store TOKUYAMA IONOTITE F at temperatures between 0 to 25°C (32 to 77°F).
- AVOID direct exposure of TOKUYAMA IONOTITE F to light and heat.

- DO NOT use Powder and Liquid after the indicated expiration date on the bottle/package.

■ **DISPOSAL**

- To dispose of TOKUYAMA IONOTITE F safely, mix the Powder and Liquid and allow the mix to harden before disposal.
- The remaining Liquid should be absorbed by an inert absorbent material such as gauze or cotton, and disposed of in accordance with local regulations.

■ **CLINICAL PROCEDURE**

- Pulp Protection
 - Glass ionomer lining material or calcium hydroxide should be applied if the cavity is in close proximity to the pulp prior to taking an impression.
 - DO NOT USE EUGENOL BASED MATERIALS for pulp protection and temporary sealing as these materials will inhibit setting/adhesion of TOKUYAMA IONOTITE F.
- Pretreatment of Prosthesis
 - Sandblast the adhesive surface thoroughly with 30- to 50-micron aluminum beads, clean the prosthesis ultrasonically, and dry. Protect the area of the metal surface not involved in cementation with a separator such as paraffin wax or petroleum jelly (such as Vaseline) if necessary.
 - Do not allow contamination of the adhesive surface. If the adhesive surface is contaminated by saliva or blood, remove contaminants thoroughly by ultrasonic cleaning.
 - If unable to sandblast the prosthesis, roughen the surface with a bur or diamond point.
 - Performing tin plating or applying METALTITE (sold separately) to the adhesive surface of the precious alloy portion of the prosthesis will improve adhesive strength and durability.
- Cleaning of Cavity/Abutment
 - Remove temporary restoration, residual temporary cement and hard deposits from the tooth by ultrasonic scalers or other cleaning instruments. Thoroughly clean all contaminants from the tooth surface to be bonded. Rinse with water and dry in the usual manner.
 - DO NOT desiccate the vital tooth. Desiccation can lead to post-operative sensitivity.
 - Root canal must be dried with absorbent points before cementing because excess water in the canal may cause premature detachment.
 - Elements listed below, which inhibit setting, should be removed from the tooth surface by thoroughly cleaning the tooth surface with alcohol, citric acid, or the application of phosphoric acid for 2 to 3 seconds before cementing.
 - Silicone oil from fit check material
 - Oil mist from handpiece
 - Saliva, blood and exudates
- Dispensing
 - Lightly shake the Powder container before dispensing. In case of one crown, dispense one level regular scoop of Powder and two drops of Liquid onto the mixing pad. In case of one inlay, dispense one level small scoop of Powder and one drop of Liquid.
- Mixing
 - Complete mixing with a plastic spatula within 30 seconds.
 - Do not use a metal spatula, because it may cause color change in the Cement due to the contamination by metal shavings.
 - Do not mix the Powder or Liquid with other brands of powder or liquid.
 - Do not return the dispensed Powder and Liquid to the bottles even if they are excessive.
- Application
 - Apply the resulting Cement to the adhesive surface of the prosthesis and tooth.
 - To keep the film thickness thin, seat the prosthesis as soon as possible after the application of cement to the tooth surface because tooth temperature accelerates the setting time. Special attention is required when seating a long-span bridge, a cast post or multiple restorations at one sitting.
- Cementing
 - Seat the prosthesis on the tooth. Waxy excess cement can be easily removed within 1 to 3 minutes after seating. Maintain normal occlusal pressure onto the prosthesis with a cotton roll for 6 or more minutes from initial placement to complete set.
 - Insufficient and delayed hardening will occur on the mixing pad due to oxygen inhibition when a thin sheet of Cement is spread across the pad. To test the progress of hardening place the Cement between a pair of slide glasses clasping with both palms of your hand or immersing the assembly in warm water. These procedures should raise the temperature close to 37°C (99°F) resulting in proper and timely set of the Cement.
- Final Adjustment
 - Follow standard procedures for occlusal check, adjustments, and final polishing.
 - DO NOT prepare a cast post abutment immediately after setting because insufficient adhesion may cause premature detachment of the post. In order to achieve sufficient adhesion to dentin, maintain appropriate pressure onto the cast post for 6 minutes and keep for 14 minutes without pressure (sum total 20 minutes is required).

- Dispensing
 - Lightly shake the Powder container before dispensing. In case of one crown, dispense one level regular scoop of Powder and two drops of Liquid onto the mixing pad. In case of one inlay, dispense one level small scoop of Powder and one drop of Liquid.
- Mixing
 - Complete mixing with a plastic spatula within 30 seconds.
 - Do not use a metal spatula, because it may cause color change in the Cement due to the contamination by metal shavings.
 - Do not mix the Powder or Liquid with other brands of powder or liquid.
 - Do not return the dispensed Powder and Liquid to the bottles even if they are excessive.
- Application
 - Apply the resulting Cement to the adhesive surface of the prosthesis and tooth.
 - To keep the film thickness thin, seat the prosthesis as soon as possible after the application of cement to the tooth surface because tooth temperature accelerates the setting time. Special attention is required when seating a long-span bridge, a cast post or multiple restorations at one sitting.
- Cementing
 - Seat the prosthesis on the tooth. Waxy excess cement can be easily removed within 1 to 3 minutes after seating. Maintain normal occlusal pressure onto the prosthesis with a cotton roll for 6 or more minutes from initial placement to complete set.
 - Insufficient and delayed hardening will occur on the mixing pad due to oxygen inhibition when a thin sheet of Cement is spread across the pad. To test the progress of hardening place the Cement between a pair of slide glasses clasping with both palms of your hand or immersing the assembly in warm water. These procedures should raise the temperature close to 37°C (99°F) resulting in proper and timely set of the Cement.
- Final Adjustment
 - Follow standard procedures for occlusal check, adjustments, and final polishing.
 - DO NOT prepare a cast post abutment immediately after setting because insufficient adhesion may cause premature detachment of the post. In order to achieve sufficient adhesion to dentin, maintain appropriate pressure onto the cast post for 6 minutes and keep for 14 minutes without pressure (sum total 20 minutes is required).

IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the dental professional to ensure the product is suitable for application before use. Specifications are subjected to change without notification.

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