Resin-based Dental Restorative Material

PALFIQUE BULK FLOW

ENGLISH
Read all information, precautions and notes before using.

■PRODUCT DESCRIPTION AND GENERAL INFORMATION

- PALFIQUE BULK FLOW is a low viscosity, light-cured, radiopaque composite resin. This low stress flowable material can be placed 4mm increments.
- PALFIQUE BULK FLOW incorporates Radical-Amplified Photopolymerization initiator technology (RAP technology). RAP technology facilitates a shortened light curing time and ample working time.
- PALFIQUE BULK FLOW contains 70% by weight (56% by volume) of spherical silica–zirconia filler and composite filler. The composite filler offers reduced polymerization shrinkage stress due to its round shape. The spherical filler contained in PALFIQUE BULK FLOW (mean particle size: 200nm, particle size range: 100 to 300nm) facilitates excellent gloss retention and wear resistance. PALFIQUE BULK FLOW contains Bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), Triethylene glycol dimethacrylate (TEGDMA), Bisphenol A polyethoxy methacrylate (Bis-MPEPP), Mequinoxil, Dibutyl hydroxyl toluene and UV absorber.
- PALFIQUE BULK FLOW is provided in SYRINGE.

■SHADE

U (Universal), A1, A2, A3, B1

- Please check with your local distributor for available shades.

■INDICATIONS

- Direct anterior and posterior restorations
- Cavity lining
- Blocking out cavity undercuts before fabrication indirect restorations
- Repair of porcelain/composite

■CONTRAINDICATIONS

PALFIQUE BULK FLOW contains methacrylic monomers. DO NOT use PALFIQUE BULK FLOW for patients allergic to or hypersensitive to methacrylic and related monomers or to any of the other ingredients.

■PRECAUTIONS

- DO NOT use PALFIQUE BULK FLOW for any purposes other than those listed in these instructions. Use PALFIQUE BULK FLOW only as directed herein.
- PALFIQUE BULK FLOW is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor use by non-dental care professionals.
- DO NOT use PALFIQUE BULK FLOW if the safety seals are broken or appear to have been tampered with.
- If PALFIQUE BULK FLOW causes an allergic reaction or hypersensitivity, discontinue its use immediately.
- Use examination gloves (plastic, vinyl, or latex) at all times when handling PALFIQUE BULK FLOW to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If PALFIQUE BULK FLOW comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- Avoid contact of PALFIQUE BULK FLOW with eyes, mucosal membrane, skin and clothing.
 - If PALFIQUE BULK FLOW comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If PALFIQUE BULK FLOW comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.
 - If PALFIQUE BULK FLOW comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
- PALFIQUE BULK FLOW should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- To avoid the unintentional ingestion of PALFIQUE BULK FLOW, do not leave it unsupervised within the reach of patients and children.
- Clean the placement instruments and brushes with alcohol after use.
- When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
- PALFIQUE BULK FLOW is available in semi-translucent shades. If applied in cavity Class III or IV direct restorations, the cured PALFIQUE BULK FLOW may appear dark. Avoid using PALFIQUE BULK FLOW for the direct restoration of Class III or IV when opacity is paramount.
- Light-cure PALFIQUE BULK FLOW extra orally to conform the time needed for complete hardening of PALFIQUE BULK FLOW with your light-curing unit prior to doing the restoration.
- This product is designed to be used at room temperature (18 - 30° C / 62 - 84° F). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

■CURING TIME AND DEPTH

PALFIQUE BULK FLOW enables a 4mm depth of cure (U shade) by light-curing for 10 seconds with Halogen or LED light-curing unit (Light intensity: 800mW/cm² or more, wavelength: 400-500nm). The depth of cure was determined based on test results performed in accordance with

section 7.10 “depth of cure” of ISO 4049:2009.

■STORAGE

- Store PALFIQUE BULK FLOW at a temperature between 0 - 25° C (32 - 77° F).
- AVOID direct exposure to light and heat.
- DO NOT use PALFIQUE BULK FLOW after the indicated date of expiration on the syringe / package.

■DISPOSAL

To safely dispose of excess PALFIQUE BULK FLOW, extrude unused portion from the syringe and light-cure before disposal.

■CLINICAL PROCEDURE

1. Cleaning

Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

2. Shade Selection

Prior to isolation of tooth, select the appropriate shade of PALFIQUE BULK FLOW using shade card provided for PALFIQUE BULK FLOW.

3. Isolation

A rubber dam is the preferred method of isolation.

4. Cavity Preparation

Prepare the cavity and rinse with water.

5. Pulp Protection

Calcium hydroxide should be applied if the cavity is in close proximity to the pulp. DO NOT USE EUGENOL-BASED MATERIALS to protect the pulp as these materials will inhibit curing PALFIQUE BULK FLOW. PALFIQUE BULK FLOW is not indicated for direct pulp capping.


6. Bonding system

Apply a light-cured bonding system according to its manufacturer’s instructions.

- Do not use self-cured bonding systems. We advise that light-cured bonding systems should be used with PALFIQUE BULK FLOW since other types of bonding systems such as self-cured or dual-cured bonding systems may not provide effective adhesion when used with light-cured composite resins. If you want to use self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.
- In case of porcelain/composite repairs be sure to condition the surface of porcelain/composite with a silane coupling agent BEFORE performing the bonding procedure.

7. Dispensing

- Hold the nozzle and turn the cap counterclockwise for removal. [1]
- Hold the nozzle and place the Dispensing Tip onto it. Turn the Dispensing Tip clockwise until it is securely locked. [2]

- Accessories and Dispensing Tip are exclusively designed for PALFIQUE BULK FLOW. Avoid using Dispensing Tips from other brands with PALFIQUE BULK FLOW.
- Dispensing Tips are also sold separately.
- Use gloves when attaching the Dispensing Tip.

8. Filling and contouring

Cavity can be directly filled from the Dispensing Tip. Increments should not exceed 4mm at a time.

- DO NOT apply excessive force to extrude the material from the syringe immediately after removal from the refrigerator.
- DO NOT mix PALFIQUE BULK FLOW with other brands of resin composite or other shades of PALFIQUE BULK FLOW to avoid incomplete cure or entrapment of air bubbles.
- In case of syringe, after extruding the paste, remove and discard the Dispensing Tip, wipe the nozzle with gauze and replace the cap immediately.

9. Curing

Light-cure each increment for at least the indicated time (please see the curing recommendations below). Light-cure each layer (up to 4mm increment) using Halogen light-curing or LED light-curing unit with a wavelength of 400-500 nm. Be sure the light-curing unit with an intensity of 800mW/cm² or more. Hold the light guide tip as close to the restorative as possible during light exposure.

CURING RECOMMENDATIONS

Shade	Increment Depth	Light Output	Curing Time
All shades	4mm	≥800mW/cm²	10 sec.

- If other brands of composite resins are layered over the cured composite, follow their instructions.

10. Finishing

Shape and polish the restoration. For finishing, use finishing discs and/ or fine finishing diamond points. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tools. For final polishing, use felt discs or cotton wheels with polishing paste, or suitable polishing tools.

IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the dental professional to ensure the product is suitable for application before use.

Specifications are subject to change without notice. When the product specification changes, the instructions and precautions may change also.

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Tokuyama Dental Corporation