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North America
Tokuyama Dental America Inc.
Tel: (877)378 3548 (Toll-Free)

ENGLISH


SOFRELINER
Read all information, precautions, notes and materials before using the product.

PRODUCT DESCRIPTION AND GENERAL INFORMATION

SOFRELINER is an addition-cured silicone chairside resilient lining material for dentures. Sofreliner has no odor or taste. It is provided in a new automic cartridge delivery system, together with an adhesive prepreg and accessories. The material cures in 5 minutes intraorally or in 20 minutes at room temperature (25°C / 77°F).

CONTENTS

- | | |
|--|-----|
| 1. One cartridge of paste containing Base (27g) and Catalyst (27g), to be used with the TOKUYAMA DISPENSER II (Sold separately)..... | x1 |
| 2. Primer (10 mL)..... | x2 |
| 3. Brushes for Primer..... | x2 |
| 4. Plastic cup for Primer..... | x2 |
| 5. Spatula for spreading SOFRELINER after it has been applied to the denture..... | x1 |
| 6. TD II Mixing tip (S) or (XS) for use with paste cartridge and the TOKUYAMA DISPENSER II..... | x10 |
| 7. Steel Bur T 101..... | x1 |
| 8. Finishing Wheel..... | x1 |
| * Note: Only the TOKUYAMA DISPENSER II can be used with SOFRELINER. | |

Refill

- | | |
|---|----|
| 1. One cartridge of paste containing Base (27g) and Catalyst (27g)..... | x1 |
| 2. Primer (10 mL)..... | x2 |
| 3. TD II Mixing tip S (50 tips)..... | x2 |
| 4. TD II Mixing tip XS (50 tips)..... | x2 |
| 5. SOFRELINER FINISHING SYSTEM..... | x1 |
| - Finishing Wheel..... | x1 |
| - Mandrel T..... | x1 |

CLINICAL APPLICATIONS FOR SOFRELINER

- Use SOFRELINER on patients with alveolar bone loss caused by sharp ridges of the alveolar bone and atrophy of the mucosa.
- Poor retention caused by extreme ridge atrophy
- Major undercuts or maxillary/mandibular tori
- Retention of overdentures before the final attachment is placed

STORAGE AND DISPOSAL

Keep SOFRELINER in a cool, dark place (0-25°C / 32-77°F) before using. After using SOFRELINER, the mixing tip must be replaced with the cartridge cap, which is attached to the cartridge at the time of purchase. DO NOT use the material if the expiration date has passed. To dispose the paste, dispense unused paste out from the cartridge and harden it before disposal. To dispose the Primer, remove the liquid from the bottle and soak the liquid in tissue paper before disposal.

CLINICAL PRECAUTIONS

- SOFRELINER is intended for use only by licensed dental care professionals.
- Use SOFRELINER only as instructed herein. DO NOT use SOFRELINER for any other application.
- DO NOT use this product if the safety seals are broken or if the contents have been tampered with.
- Primer contains methylene chloride, which is highly volatile. Avoid long-term exposure, ventilate the room during use, and close the cap tightly after use.
- SOFRELINER is a silicone-based product. DO NOT use SOFRELINER on patients with a history of sensitivity to silicone. Always determine whether your patient is allergic to, or hypersensitive to silicone-based products. If a patient is uncertain whether or not he/she is sensitive to silicone based products, such as SOFRELINER, always perform a sensitivity test on the patient before using this product. Allow a small amount of SOFRELINER paste to set on a dispensing dish. Use this set material as a test patch by placing the material with a sticking plaster on the patient's cheek at least 24 hours, monitor the patient's reaction for sensitivity. Advise all patients of the possibility of sensitivity to silicone, and discuss the signs of silicone sensitivity with them. If a patient shows any signs of sensitivity, patient should be instructed to immediately contact their dentist for removal of SOFRELINER from the denture.
- The use of examination gloves (plastic, vinyl) is recommended at all times to avoid the possibility of allergic reactions from SOFRELINER. The use of latex gloves may inhibit the cure of SOFRELINER. If SOFRELINER Primer comes in contact with the examination gloves, immediately dispose of the examination gloves and wash hands thoroughly with water.
- AVOID CONTACT WITH EYES. In case of contact with eyes, thoroughly flush eyes with water, and immediately contact a physician.
- Avoid contact of uncured chemical with skin. In case of contact with skin, clean effected area with rubbing alcohol, then thoroughly flush with water.
- Avoid contact of unpolymerized materials with clothing. In case of contact with clothing, do not wipe the material. Allow it to set and peel it off.
- DO NOT allow excess material to flow down the patient's trachea or esophagus. Aspirated or ingested SOFRELINER paste can cause asphyxiation or other injury, even though SOFRELINER paste is non-toxic. Immediately remove any excess lining material from the palatal area to avoid flow of material down the patient's throat. SOFRELINER is not radioopaque, and X-rays will not detect aspirated or ingested material.
- DO NOT allow excess material to flow down the patient's trachea or esophagus. Aspirated or ingested SOFRELINER paste can cause asphyxiation or other injury, even though SOFRELINER paste is non-toxic. Immediately remove any excess lining material from the palatal area to avoid flow of material down the patient's throat. SOFRELINER is not radioopaque, and X-rays will not detect aspirated or ingested material.

- Apply SOFRELINER Primer to the bond surface of the acrylic. Use only the Primer supplied with SOFRELINER, and use only a thin coat of Primer. DO NOT apply SOFRELINER Primer to the patient's mouth or gums. SOFRELINER Primer should be applied only to the denture and always with the caution outside of the patient's mouth.
- If the SOFRELINER Primer has thickened, do not use it.
- Always use a mixing tip to dispense.
- For a thick surface, allow the SOFRELINER Paste dispensed on to the denture to stiffen slightly before inserting the denture into the patient's mouth.
- AVOID CONTACT WITH EYES. In case of contact with eyes, thoroughly clean the impression material from the denture surface before applying SOFRELINER. This is particularly important when using the laboratory lining method with a lining jig.
- Soaking the denture in heavy bleach cleaning solutions may alter the color of the SOFRELINER, but it can be used continuously.
- Instruct the patient to clean the denture as described in the following section.

PATIENT PRECAUTIONS

- Patients should be informed of the following to enhance the longevity, durability and safety of a denture lined with SOFRELINER material:
- SOFRELINER material should not be scrubbed with a denture brush or toothbrush. The best way to clean SOFRELINER is with cotton or a soft piece of gauze. While running cold water on the SOFRELINER surface, wipe the surface with the cotton or gauze dipped in soap. Use a gentle wiping motion to remove all debris and food deposits from the SOFRELINER material. The tooth surface and hard acrylic denture base may be cleaned with a regular denture brush by scrubbing with soap, denture toothpaste or regular denture cleanser at the dentist's discretion.
 - Soak the denture at least once a day, preferably at bedtime, in an oxygenating denture cleanser (i.e., Polident, Efferdent, etc.). The denture should not be cleaned with oxygenating cleanser within the patient's mouth. Color fading may occur with long term use, depending upon the patient's denture maintenance.
 - Soaking the denture in heavy bleach cleaning solutions may alter the color of the SOFRELINER, but it can be used continuously.
 - The denture should be stored in water when not being worn.
 - The lined denture should be inspected and ultrasonically cleaned by a dentist every three months, or when plaque deposit is noted.
 - If plaque is routinely present on the denture, the denture should be soaked at least every other day, in plaque control mouthwash.
 - The SOFRELINER must be inspected frequently for any damage. If the liner is disturbed, broken, cracked, loose or damaged in any way, the patient should inform his/her dentist immediately. A denture with disturbed, broken, cracked, loose or damaged liner should not be worn. Use of a damaged or broken denture could result in choking, aspiration, asphyxiation or other injury.
 - In the event of any skin sensitivity to silicone, use of the denture should be discontinued, and the patient should contact his/her dentist immediately.

- The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure the product is suitable for the appropriate application and use.
- Specifications are subjected to change without notification.

CLINICAL APPLICATION OF SOFRELINER INTRAORAL LINING METHOD

- To achieve the maximum effect using SOFRELINER, the optimum thickness of the lining material should be 2-3 millimeters. The set of the material will not be affected if SOFRELINER is applied thicker than 3 millimeters in some areas. Areas which will require lining thicker than 3 millimeters should first be lined with TOKUYAMA REBASE II, Tokuyama's chairside hard denture lining material (sold separately), to obtain the proper thickness for SOFRELINER.
- Remove the denture from the patient's mouth and clean the area to be lined with soap and water. Then air-dry the surface of the denture thoroughly.
- Grind the tissue side of the denture to the required thickness with a bur. Bevel the borders of the denture to allow at least 2 millimeters of material to roll on the borders. Avoid setting the marginal line on the denture border.
- After blowing the denture with strong air stream, roll the Primer into the plastic cup included in the SOFRELINER kit. Using the brush included in the SOFRELINER kit, apply the Primer evenly in a thin coat on the surface of the denture where SOFRELINER is to be lined. Do not miss any spot where the Primer should be applied, since this may cause poor adhesion.
- Note: Attach the Drip-cap included in the SOFRELINER kit onto the rim of the Primer bottle to facilitate pouring.
- Load the SOFRELINER paste cartridge into a TOKUYAMA DISPENSER II (sold separately).
- Evenly dispense SOFRELINER paste on the tissue side of the denture. Spread the paste evenly with a spatula across all areas where the Primer was applied.
- Discontinue any further SOFRELINER paste on the patient's mouth. In 1/2 minutes after dispensing SOFRELINER onto the denture. Place the denture in the patient's mouth, and have the patient occlude very gently into centric occlusion to the proper vertical dimension. Check for any excess SOFRELINER in the patient's mouth, and remove immediately if any is present. DO NOT allow excess material to flow down the patient's trachea or esophagus. Aspirated or ingested SOFRELINER paste can cause asphyxiation or other injury, even though SOFRELINER paste is non-toxic. If the material flows towards the patient's trachea, immediately wipe excess material from the patient's mouth with a

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North America
Tokuyama Dental America Inc.
Tel: (877)378 3548 (Toll-Free)

- cotton ball or gauze. If, at any time, the patient feels SOFRELINER material or any substance flowing toward his/her throat, they should immediately alert the dentist.
Note: To line upper dentures with SOFRELINER, the denture from the posterior side should be removed. After each use, the patient should resume occluding gently.
- Put the checks very gently to border mold the vestibule area while the patient is gently occluding.
 - Allow the denture to completely set in the patient's mouth while continuing to monitor the patient for the presence of excess SOFRELINER. If patient salivates, suction the excess SOFRELINER material on the denture.
 - SOFRELINER should be allowed to set intraorally for at least 5 minutes.
 - Remove the denture from the patient's mouth, and carefully trim any excess SOFRELINER material from the border areas with a knife or a sharp scalpel. DO NOT over trim the area where the bevel was made on the denture. DO NOT leave any loose SOFRELINER material on the denture.
 - *Note: An additional layer of SOFRELINER can be added to a newly lined layer of SOFRELINER without the use of Primer only if it is performed immediately after the initial lining when the denture has been washed carefully and dried completely. SOFRELINER Paste will not bond to completely cured SOFRELINER.

FINISHING

- Use the Steel Bur T 101 included in the SOFRELINER Kit. When necessary, use the Finishing Wheel included in the SOFRELINER kit for the final finishing.
- OTHER IMPORTANT CONSIDERATIONS FOR THE SUCCESSFUL USE OF SOFRELINER
In order to enhance bonding on old or worn dentures, line the denture with TOKUYAMA REBASE II prior to lining with SOFRELINER. It is important to be sure that the tissue side of the denture is ground deeply enough, so that the relining of SOFRELINER will not be too high on the patients' occlusion after the lining with TOKUYAMA REBASE II.

CLINICAL APPLICATION OF SOFRELINER LABORATORY LINING METHOD

- EXAMINATION OF THE DENTURE
If the denture is worn or old, line the denture with TOKUYAMA REBASE II. The TOKUYAMA REBASE II material should be applied in the appropriate thickness so that when the new material is ground down, a newly lined layer of TOKUYAMA REBASE II is left on the denture bearing area allowing SOFRELINER to adhere to it. In this case, make certain that the tissue side of the denture is ground deeply enough, so that the relined SOFRELINER will not be too high on the patients' occlusion after the lining with TOKUYAMA REBASE II.
- FABRICATION OF A GYPSUM MODEL
Take an impression with the patient biting on the denture, then pour a stone model on the impression.
 - Make a stone index of the teeth and mount the denture on the index of a lining jig.
 - Prepare the denture for the thickness of the lining material by grinding the lining area with an acrylic bur.
 - Replace the denture on the index and position the denture on the lining jig.
 - Apply a separating medium on the ridge area of the gypsum model so the lining material will not adhere to the gypsum model.
- Apply Primer to the entire area of the denture to be lined.
- Place the cartridge in the dispenser and apply the paste onto the surface where the Primer has been applied.
- Place the denture onto its gypsum stone model, then press it with the lining jig.
- Allow denture to sit for 20 minutes at room temperature (23°C / 73.4°F) on the lining jig. *1
- After SOFRELINER has been cured, remove the denture from the gypsum model, and carefully trim any excess SOFRELINER material from the border areas with a knife or a sharp scalpel. DO NOT over trim the area where the bevel was made on the denture. DO NOT leave any loose SOFRELINER material on the denture. *2
- *1 Always maintain the denture in the lining jig for 20 minutes or more at room temperature (23°C / 73.4°F) in order to obtain maximum room temperature set.
- *2 An additional layer of SOFRELINER can be added to a newly lined layer of SOFRELINER without the use of Primer, only if it is performed immediately after the initial lining. SOFRELINER paste will not bond to completely cured SOFRELINER.

FINISHING

- Use the Steel Bur T 101 included in the SOFRELINER kit for the final finishing.
- ### HELPFUL HINTS FOR USING SOFRELINER
- If there is a great amount of ridge resorption, line the old or worn denture with a new acrylic plaster on the patient's cheek to enhance the bond between the denture and SOFRELINER.
 - Use SOFRELINER Primer to bond the surface of the acrylic. Use only the Primer supplied with SOFRELINER, and use only a thin coat of Primer. DO NOT apply SOFRELINER Primer to the patient's mouth or gums. SOFRELINER Primer should be applied only to the denture and always with the caution outside of the patient's mouth.
 - If the SOFRELINER Primer has thickened, do not use it.
 - Always use a mixing tip to dispense.
 - For a thick surface, allow the SOFRELINER Paste dispensed on to the denture to stiffen slightly before inserting the denture into the patient's mouth.
 - AVOID CONTACT WITH EYES. In case of contact with eyes, thoroughly clean the impression material from the denture surface before applying SOFRELINER. This is particularly important when using the laboratory lining method with a lining jig.
 - Soaking the denture in heavy bleach cleaning solutions may alter the color of the SOFRELINER, but it can be used continuously.
 - Instruct the patient to clean the denture as described in the following section.

PATIENT PRECAUTIONS

- Patients should be informed of the following to enhance the longevity, durability and safety of a denture lined with SOFRELINER material:
- SOFRELINER material should not be scrubbed with a denture brush or toothbrush. The best way to clean SOFRELINER is with cotton or a soft piece of gauze. While running cold water on the SOFRELINER surface, wipe the surface with the cotton or gauze dipped in soap. Use a gentle wiping motion to remove all debris and food deposits from the SOFRELINER material. The tooth surface and hard acrylic denture base may be cleaned with a regular denture brush by scrubbing with soap, denture toothpaste or regular denture cleanser at the dentist's discretion.
 - Soak the denture at least once a day, preferably at bedtime, in an oxygenating denture cleanser (i.e., Polident, Efferdent, etc.). The denture should not be cleaned with oxygenating cleanser within the patient's mouth. Color fading may occur with long term use, depending upon the patient's denture maintenance.
 - Soaking the denture in heavy bleach cleaning solutions may alter the color of the SOFRELINER, but it can be used continuously.
 - The denture should be stored in water when not being worn.
 - The lined denture should be inspected and ultrasonically cleaned by a dentist every three months, or when plaque deposit is noted.
 - If plaque is routinely present on the denture, the denture should be soaked at least every other day, in plaque control mouthwash.
 - The SOFRELINER must be inspected frequently for any damage. If the liner is disturbed, broken, cracked, loose or damaged in any way, the patient should inform his/her dentist immediately. A denture with disturbed, broken, cracked, loose or damaged liner should not be worn. Use of a damaged or broken denture could result in choking, aspiration, asphyxiation or other injury.
 - In the event of any skin sensitivity to silicone, use of the denture should be discontinued, and the patient should contact his/her dentist immediately.

The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure the product is suitable for the appropriate application and use.

Specifications are subjected to change without notification.

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
- SOFRELINER ist zur Anwendung durch Fachpersonal vorgesehen.
- SOFRELINER nur für den vorgesehenen Zweck verwenden.
- Wenden Sie das Produkt nicht mehr an, wenn die Versiegelung beschädigt oder mit dem Inhalt ungeschickmäßig umgegangen worden ist.
- EXAMINATION OF THE DENTURE
If the denture is worn or old, line the denture with TOKUYAMA REBASE II. The TOKUYAMA REBASE II material should be applied in the appropriate thickness so that when the new material is ground down, a newly lined layer of TOKUYAMA REBASE II is left on the denture bearing area allowing SOFRELINER to adhere to it. In this case, make certain that the tissue side of the denture is ground deeply enough, so that the relined SOFRELINER will not be too high on the patients' occlusion after the lining with TOKUYAMA REBASE II.
- FABRICATION OF A GYPSUM MODEL
Take an impression with the patient biting on the denture, then pour a stone model on the impression.
 - Make a stone index of the teeth and mount the denture on the index of a lining jig.
 - Prepare the denture for the thickness of the lining material by grinding the lining area with an acrylic bur.
 - Replace the denture on the index and position the denture on the lining jig.
 - Apply a separating medium on the ridge area of the gypsum model so the lining material will not adhere to the gypsum model.
- Apply Primer to the entire area of the denture to be lined.
- Place the cartridge in the dispenser and apply the paste onto the surface where the Primer has been applied.
- Place the denture onto its gypsum stone model, then press it with the lining jig.
- Allow denture to sit for 20 minutes at room temperature (23°C / 73.4°F) on the lining jig. *1
- After SOFRELINER has been cured, remove the denture from the gypsum model, and carefully trim any excess SOFRELINER material from the border areas with a knife or a sharp scalpel. DO NOT over trim the area where the bevel was made on the denture. DO NOT leave any loose SOFRELINER material on the denture. *2
- *1 Always maintain the denture in the lining jig for 20 minutes or more at room temperature (23°C / 73.4°F) in order to obtain maximum room temperature set.
- *2 An additional layer of SOFRELINER can be added to a newly lined layer of SOFRELINER without the use of Primer, only if it is performed immediately after the initial lining. SOFRELINER paste will not bond to completely cured SOFRELINER.

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