Silicone-based soft denture liner

SOFRELINER TOUGH M Medium

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ESPAÑOL
PORTUGUÊS

Carefully read all of these instructions before using SOFRELINER TOUGH M.

- PRODUCT DESCRIPTION AND GENERAL INFORMATION**
SOFRELINER TOUGH M is an addition-cured silicone chairside soft lining material for removable dentures.
- SOFRELINER TOUGH M Kit contains PASTE, PRIMER and other accessories such as mixing tips, drip nozzle, drip cap, brushes, SHAPE ADJUSTMENT POINT (White), FINISHING POINT (Brown).
- PRIMER is used solely for bonding to the acrylic surfaces. It cannot be used with nylon, metal and silicone surfaces. PASTE is dispensed in automx cartridge delivery system. The material sets in 5 minute intraorally or 20 minutes at room temperature (23° C / 73.4° F) and has no odor or taste.

INDICATIONS

- For pain caused by sharp ridges of alveolar bone and mucous atrophy.
- Pain reduction caused by extreme ridge atrophy.
- Major undercuts or maxillary / mandibular tori.
- Retention of overdentures before the final attachment is placed.

CONTRAINDICATIONS

SOFRELINER TOUGH M is contraindicated in patients allergic to or hypersensitive to silicone and platinum compounds and should NOT be used on such patients.

PRECAUTIONS

- SOFRELINER TOUGH M is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
- Use SOFRELINER TOUGH M only as instructed herein. DO NOT use SOFRELINER TOUGH M for any application other than those noted in these instructions.
- DO NOT use SOFRELINER TOUGH M if the safety seal is broken or if it appears the contents have been tampered with.

- PRIMER is flammable, which is a highly flammable and volatile solvent. Do not expose the PRIMER to open flame. Use in a ventilated room during use and close the cap tightly after use.
- If patient experiences an allergic reaction or sensitivity to SOFRELINER TOUGH M, such as a dermatitis outbreak, remove the prosthesis and discontinue its use immediately.
- Wear non-latex examination gloves at all times while using SOFRELINER TOUGH M to avoid the possibility of allergic reactions. Latex gloves may also inhibit the setting of SOFRELINER TOUGH M. If SOFRELINER TOUGH M PRIMER comes in contact with the examination gloves, immediately dispose of the examination gloves and wash hands thoroughly with water.
- Avoid contact with eyes. In case of contact with eyes, thoroughly flush eyes with water, and immediately contact a physician.
- Avoid contact of unset chemical with skin. In case of contact with skin, clean affected area with rubbing alcohol, then thoroughly flush with water.
- Avoid contact of unset materials with clothing. In case of contact with clothing, do not wipe the material. Allow it to set and then peel it off.
- DO NOT allow excess materials to flow down the patient's trachea or esophagus. Aspiration or ingested SOFRELINER TOUGH M PASTE can cause asphyxiation or other injury. Immediately remove any excess lining material from the palatal area and avoid material flowing down the patient's throat. SOFRELINER TOUGH M is not radiopaque, and radiographs will not detect aspiration or ingested material. In case of an elderly patient, radiopacity is compromised, set the chair in an upright position and instruct the patient to breathe through the nose in order to prevent the material from entering the airway.
- DO NOT mix SOFRELINER TOUGH M with any materials not noted in these instructions.

- After using SOFRELINER TOUGH M cartridge, the mixing tip must be replaced with the cartridge cap that comes with the cartridge at the time of purchase. Wipe off the excess paste on the tip of cartridge with a clean cloth or tissue each time before using a new mixing tip or replacing the cartridge cap.
- Always clean the spatula, brushes and cups with alcohol after each use to avoid cross contamination.
- Instruct patients in accordance with the PATIENT PRECAUTIONS contained in the last section of these instructions.
- To avoid inhaling dust while finishing, always wear a mask recommended by federal and/or local regulations.
- Read and follow the manufacturers' instructions for all products mentioned in this instruction manual that are not included in this kit.
- As much as possible, avoid contact SOFRELINER TOUGH M with oral mucous membrane other than the region to be lined. Instruct patients to rinse their mouth with adequate amount of water after the treatment.
- DO NOT place the PRIMER in direct contact with plastic materials, contact may dissolve some types of plastic.

- NOTES FOR USE**
 - In the case of a poor fitting denture, first line the denture with an acrylic material, such as TOKUYAMA REBASE II, chairside hard denture lining material, or any other laboratory/indirect procedures to improve it.
 - In the case of an old worn denture, the surface to be lined should be lined with new acrylic materials absorbed in water at room temperature (23° C / 73.4° F).
 - If the SOFRELINER TOUGH M PRIMER becomes too thick due to volatilization, do not use it as the adhesive strength may be reduced.

- Regarding the brush for the primer, wipe off the tip with alcohol-soaked cotton swab or paper after each use. Contaminated brush can lower the adhesion of SOFRELINER TOUGH M.
- Do not mix pastes by hand, doing so may cause air bubbles to be incorporated into the paste, which could result in deterioration, staining and unpleasant odors.
- Avoid contact SOFRELINER TOUGH M with the materials below, as they may inhibit the setting of SOFRELINER TOUGH M. Silicone-based impression materials, enamel containing materials, polysulfated rubber impression materials, natural rubber gloves, other manufacturer's soft lining materials and uncured acrylic resins.
- The material sets in 5 minutes intraorally or 20 minutes at room temperature (23° C / 73.4° F). After 24 hours the liner turns slightly firm. The finishing procedure can be performed immediately after the removal of the newly lined denture from the mouth or a gypsum model.
- Soaking the denture in solutions containing heavy bleach may alter the color of SOFRELINER TOUGH M; however, it will not affect its performance.
- Instruct patient to clean the denture as described in PATIENT PRECAUTIONS.
- Since silicone-based fit checking material sticks to SOFRELINER TOUGH M, apply Vaseline on the surface prior to use.
- If there is a fracture in the finished liner, completely remove the liner and repeat the procedure.
- Refer to following tools for shape adjustment and finishing.

Stages	Tools	Shapes
Trimming	Sharp scalpel	Slightly rounded blade
Shape Adjustment	SHAPE ADJUSTMENT POINT (White)	Round tapered end
Finishing	FINISHING POINT (Brown)	Round tapered end
Relief	DIAMOND ROUND BUR (sold separately)	Round
STORAGE		
1. Store SOFRELINER TOUGH M at a temperature between 0 – 25° C (32 – 77° F).		
2. Avoid direct exposure of SOFRELINER TOUGH M to heat, open flame, humidity and sunlight.		
3. Do not use SOFRELINER TOUGH M after the indicated date of expiration.		
4. DO NOT leave SOFRELINER TOUGH M unnecessarily within the reach of patients and children to avoid being swallowed by mistake.		
DISPOSAL		
1. To dispose of the PASTE, dispense unused paste from the cartridge and allow it to set prior to disposal.		
2. The remaining PRIMER liquid should be absorbed in an inert absorbent material such as gauze or cotton, and disposed of in accordance with local regulations.		
CLINICAL APPLICATION OF SOFRELINER TOUGH M INTRAORAL LINGING METHOD		
In most cases the optimum thickness of SOFRELINER TOUGH M should be 1-2 millimeter. If a thickness of over 2 millimeters is required, the area should first be lined with TOKUYAMA REBASE II, chairside hard denture liner, to obtain sufficient denture thickness. After the denture is lined, grind the tissue surface of the denture with a cross cut carbide bur (for acrylic resin) to the proper depth (See Figure 2).		
1. Preparation		

- Check intraoral condition by inspection and palpation. Remove denture from patient's mouth and clean the area to be lined with soap and water. Rinse thoroughly and air-dry the denture surface completely. Grind the tissue side of the denture with a cross cut carbide bur to the proper depth. Bevel the borders of the denture to allow 1 millimeter of material to roll over the borders (See Figure 1). Remove the grinding dust with a dry brush or dry oil-free air.
- Note: When the border is placed within the flange of denture base, grind the tissue side of the denture with round bur (for acrylic resin) to the proper depth (See Figure 2).


Figure 1




Figure 2

- Dispensing the PRIMER**
Open the bottle of PRIMER. Attach the drip nozzle to the bottle. Dispense an adequate amount of PRIMER into the drip-cup. Replace the cap, closing the bottle tight immediately after dispensing.

- Applying the PRIMER**
Using the brush for the PRIMER, apply a thin, even coat of PRIMER on the denture surface where SOFRELINER TOUGH M is to be applied. Do not miss any areas where the PRIMER should be applied since this may cause poor adhesion. Then air-dry the surface thoroughly until no visible liquid remains.

- Dispensing the PASTE**
Load the SOFRELINER TOUGH M PASTE cartridge into a dispenser (sold separately). Open the cartridge cap. PASTE. Attach the mixing tip to the cartridge by aligning the notches, push it onto the cartridge, rotating the mixing tip clockwise 90 degrees.

- Applying the PASTE**
Extrude the mixed paste evenly on the tissue side of the denture. Using a spatula spread the paste evenly across all surfaces where the PRIMER was applied.

Note: If any part of the primed surface is not covered with paste in this stage, the area (CROSS CUT CARBIDE BUR) should be resurfaced by the dentist.

- PRIMER is used solely for bonding to the acrylic surfaces. It cannot be used with nylon, metal and silicone surfaces. PASTE is dispensed in automx cartridge delivery system. The material sets in 5 minute intraorally or 20 minutes at room temperature (23° C / 73.4° F) and has no odor or taste.

INDICATIONS

- For pain caused by sharp ridges of alveolar bone and mucous atrophy.
- Pain reduction caused by extreme ridge atrophy.
- Major undercuts or maxillary / mandibular tori.
- Retention of overdentures before the final attachment is placed.

CONTRAINDICATIONS

SOFRELINER TOUGH M is contraindicated in patients allergic to or hypersensitive to silicone and platinum compounds and should NOT be used on such patients.

PRECAUTIONS

- SOFRELINER TOUGH M is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
- Use SOFRELINER TOUGH M only as instructed herein. DO NOT use SOFRELINER TOUGH M for any application other than those noted in these instructions.
- DO NOT use SOFRELINER TOUGH M if the safety seal is broken or if it appears the contents have been tampered with.

- PRIMER is flammable, which is a highly flammable and volatile solvent. Do not expose the PRIMER to open flame. Use in a ventilated room during use and close the cap tightly after use.
- If patient experiences an allergic reaction or sensitivity to SOFRELINER TOUGH M, such as a dermatitis outbreak, remove the prosthesis and discontinue its use immediately.
- Wear non-latex examination gloves at all times while using SOFRELINER TOUGH M to avoid the possibility of allergic reactions. Latex gloves may also inhibit the setting of SOFRELINER TOUGH M. If SOFRELINER TOUGH M PRIMER comes in contact with the examination gloves, immediately dispose of the examination gloves and wash hands thoroughly with water.
- Avoid contact with eyes. In case of contact with eyes, thoroughly flush eyes with water, and immediately contact a physician.
- Avoid contact of unset chemical with skin. In case of contact with skin, clean affected area with rubbing alcohol, then thoroughly flush with water.
- Avoid contact of unset materials with clothing. In case of contact with clothing, do not wipe the material. Allow it to set and then peel it off.
- DO NOT allow excess materials to flow down the patient's trachea or esophagus. Aspiration or ingested SOFRELINER TOUGH M PASTE can cause asphyxiation or other injury. Immediately remove any excess lining material from the palatal area and avoid material flowing down the patient's throat. SOFRELINER TOUGH M is not radiopaque, and radiographs will not detect aspiration or ingested material. In case of an elderly patient, radiopacity is compromised, set the chair in an upright position and instruct the patient to breathe through the nose in order to prevent the material from entering the airway.
- DO NOT mix SOFRELINER TOUGH M with any materials not noted in these instructions.

- After using SOFRELINER TOUGH M cartridge, the mixing tip must be replaced with the cartridge cap that comes with the cartridge at the time of purchase. Wipe off the excess paste on the tip of cartridge with a clean cloth or tissue each time before using a new mixing tip or replacing the cartridge cap.
- Always clean the spatula, brushes and cups with alcohol after each use to avoid cross contamination.
- Instruct patients in accordance with the PATIENT PRECAUTIONS contained in the last section of these instructions.
- To avoid inhaling dust while finishing, always wear a mask recommended by federal and/or local regulations.
- Read and follow the manufacturers' instructions for all products mentioned in this instruction manual that are not included in this kit.
- As much as possible, avoid contact SOFRELINER TOUGH M with oral mucous membrane other than the region to be lined. Instruct patients to rinse their mouth with adequate amount of water after the treatment.
- DO NOT place the PRIMER in direct contact with plastic materials, contact may dissolve some types of plastic.

- NOTES FOR USE**
 - In the case of a poor fitting denture, first line the denture with an acrylic material, such as TOKUYAMA REBASE II, chairside hard denture lining material, or any other laboratory/indirect procedures to improve it.
 - In the case of an old worn denture, the surface to be lined should be lined with new acrylic materials absorbed in water at room temperature (23° C / 73.4° F).
 - If the SOFRELINER TOUGH M PRIMER becomes too thick due to volatilization, do not use it as the adhesive strength may be reduced.

- Regarding the brush for the primer, wipe off the tip with alcohol-soaked cotton swab or paper after each use. Contaminated brush can lower the adhesion of SOFRELINER TOUGH M.
- Do not mix pastes by hand, doing so may cause air bubbles to be incorporated into the paste, which could result in deterioration, staining and unpleasant odors.
- Avoid contact SOFRELINER TOUGH M with the materials below, as they may inhibit the setting of SOFRELINER TOUGH M. Silicone-based impression materials, enamel containing materials, polysulfated rubber impression materials, natural rubber gloves, other manufacturer's soft lining materials and uncured acrylic resins.
- The material sets in 5 minutes intraorally or 20 minutes at room temperature (23° C / 73.4° F). After 24 hours the liner turns slightly firm. The finishing procedure can be performed immediately after the removal of the newly lined denture from the mouth or a gypsum model.
- Soaking the denture in solutions containing heavy bleach may alter the color of SOFRELINER TOUGH M; however, it will not affect its performance.
- Instruct patient to clean the denture as described in PATIENT PRECAUTIONS.
- Since silicone-based fit checking material sticks to SOFRELINER TOUGH M, apply Vaseline on the surface prior to use.
- If there is a fracture in the finished liner, completely remove the liner and repeat the procedure.
- Refer to following tools for shape adjustment and finishing.

- 1. Removal**
Remove a large portion of SOFRELINER TOUGH M material by cutting it off with a sharp instrument. Remove the remaining material by grinding with a diamond round bur. Note: SILICONE REMOVER (sold separately) easily removes old silicone denture liner from the acrylic denture base.
- CLINICAL APPLICATION OF SOFRELINER TOUGH M LABORATORY LINING METHOD**
1. Examination of mouth
Check intraoral condition by inspection and palpation. Examine the area where needs relief.
- 2. Examination of the denture**
If the denture is worn or old, line the denture first using TOKUYAMA REBASE II, chairside hard denture lining material. The hard lining material should be applied in the appropriate thickness so that when the new material is ground down according to the procedure in step 4, a newly lined layer of hard lining material is left on the denture bearing area, allowing SOFRELINER TOUGH M to adhere to it.
- 3. Fabrication of a wax model**
Take an impression with the patient biting on the denture, then pour a stone slurry onto the impression. Make a stone index of the teeth and mount the denture on the index of a lining jig. Relieve the stone model with wax in the region where may cause pain in the mouth. Apply a separating medium on the ridge area of the gypsum model so the lining material will adhere to it.
- 4. Preparation**
Grind the tissue side of the denture with a cross cut carbide bur to the proper depth. Bevel the borders of the denture to allow 1 millimeter of material to roll over the borders (See Figure 1). Remove the grinding dust with a dry brush or dry oil-free air.
- Note: Use of a round bur makes it easy to remove a uniform depth of resin. The diameter of the bur can be used as a guide for depth.

- Dispensing the PRIMER**
Open the bottle of PRIMER. Attach the drip nozzle to the bottle. Dispense an adequate amount of PRIMER into the drip-cup. Replace the cap, closing the bottle tight immediately after dispensing.
- Applying the PRIMER**
Using the brush for the PRIMER, apply a thin, even coat of PRIMER on the denture surface where SOFRELINER TOUGH M is to be lined. Do not miss any areas where the PRIMER should be applied since this may cause poor adhesion. Then air-dry the surface thoroughly until no visible liquid remains.

- Dispensing the PASTE**
Load the SOFRELINER TOUGH M PASTE cartridge into a dispenser (sold separately). Open the cartridge cap. PASTE. Attach the mixing tip to the cartridge by aligning the notches, push it onto the cartridge rotating the mixing tip clockwise 90 degrees.

- Applying the PASTE**
Extrude the mixed paste onto the surface where the PRIMER has been applied. Position the paste both on the mucosal side of denture and stone model.

- Trimming**
Reattach the denture onto its gypsum stone model and press it with the lining jig. Leave undisturbed for 20 minutes at room temperature (23° C / 73.4° F).

- Finishing**
Note: If the set material is sticky, soak the denture in warm water at 40° C to 50° C (104° F to 122° F) for 5 minutes with the lining jig.

- Relief**
Remove the denture from the gypsum model. Carefully trim all excess SOFRELINER TOUGH M material from the border areas with a sharp scalpel (See Figure 3).

- Finishing**
For shaping, use the SHAPE ADJUSTMENT POINT (White), by applying light pressure not to exceed 15,000 rpm allowing the point to move in a small circle. For finishing, use the FINISHING POINT (Brown), by applying light pressure not to exceed 4,000 rpm.

- Note: Rotate both the SHAPE ADJUSTMENT POINT and the FINISHING POINT from the lined area to the denture so as not to break the margin region (See Figure 4).

- Use DIAMOND ROUND BUR (sold separately) or any suitable tool for relief.

PATIENT PRECAUTIONS

Patients should be informed of the following instructions to enhance the longevity, durability and safety of a denture lined with SOFRELINER TOUGH M:

- SOFRELINER TOUGH M lining material should be scrubbed with a denture brush or toothbrush. The best way to clean SOFRELINER TOUGH M is with a sponge or a soft piece of gauze. While running cold water on the SOFRELINER TOUGH M surface, wipe the surface with sponge or gauze saturated with soap. Use a gentle wiping motion to remove all debris and food deposits from the SOFRELINER TOUGH M lining material. The tooth surface and hard acrylic denture base may be cleaned with a regular denture brush by scrubbing with soap, denture toothpaste or regular denture cleansers at the dentist's discretion.

- Soak the denture at least once a day, preferably at bedtime, in an oxygenating denture cleaner. The denture should not be cleaned with an oxygenating cleaner when the patient's mouth. Color fading may occur with long term tissue color depending on the patient's denture maintenance. Use a denture cleaner according to its manufacturer instructions.

- Soaking the denture in heavy bleach cleaning solutions may alter the color of the SOFRELINER TOUGH M, but it will not affect its performance.

- The denture should be stored in water when not being worn.

- The lined denture should be inspected ultrasonically cleaned by a dentist every three months, or whenever plaque deposit is noted.

- The SOFRELINER TOUGH M should be inspected frequently for any damage. A denture with disturbed, broken, cracked, loose or damaged liner should not be worn.

- Use of a damaged or broken denture could result in choking, aspiration, asphyxiation or other injury.

- In the event of any sign of sensitivity to silicone, use of the denture should be discontinued, and the patient should contact his/her dentist immediately.

IMPORTANT NOTE: The manufacturer of SOFRELINER TOUGH M is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure that this product is suitable for an appropriate application before use.

Product specifications of SOFRELINER TOUGH M are subjected to change without notification. When the product specifications change, the instruction and precaution may change also.

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