Resin-based Dental Restorative Material

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QUICK-CURE

ENGLISH

Read all information, precautions and notes before using.

■ PRODUCT DESCRIPTION AND GENERAL INFORMATION

- 1) ESTELITE POSTERIOR is a light-cured, radiopaque composite resin for use in posterior restorations. ESTELITE POSTERIOR contains 84% by weight (70% by volume) of silica-zirconia filler (mean particle size: 2 μm , particle size range: 0.1 to 10 μm). A high filler load offers decreased polymerization shrinkage. The monomer matrix contains Bis-GMA, Triethylene glycol dimethacrylate (TEGMDA) and Bisphenol A polyethoxy methacrylate (Bis-MPEPP).
- 2) ESTELITE POSTERIOR implements Radical-Amplified Photopolymerization initiator technology (RAP technology). Curing time for ESTELITE POSTERIOR can be shortened by one-third using a light curing unit with a wavelength range of 400 to 500nm (peak: 470nm) compared to our conventional resin composite. Please see the table depicting the relationship between curing time and increment depth (INDICATIONS FOR FILLING AND CURING).
- 3) RAP technology facilitates a longer working time for ESTELITE POSTERIOR.
- 4) ESTELITE POSTERIOR is provided in SYRINGE.

■ SHADE

- PCE, PA1, PA2, PA3.
- Every shade provided for ESTELITE POSTERIOR has adequate translucency designed for posterior restoration.
- PCE serves as a top layer for the multi-shade layering technique.
- Tokuyama Dental Corporation does not provide ESTELITE POSTERIOR Shade Guide. Please make a mockup shade guide by using ESTELITE POSTERIOR.

■ INDICATIONS

- Direct posterior restorations including the occlusal surface
- Repair of porcelain/composite

■ CONTRAINDICATIONS

ESTELITE POSTERIOR contains methacrylic monomers. DO NOT use ESTELITE POSTERIOR for patients allergic to or hypersensitive to methacrylic and related monomers.

■ PRECAUTIONS

- 1) DO NOT use ESTELITE POSTERIOR for any purpose other than those listed in these instructions. Use ESTELITE POSTERIOR only as directed herein.
2) ESTELITE POSTERIOR is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
3) DO NOT use ESTELITE POSTERIOR if the safety seals are broken or appear to have been tampered with.
4) If ESTELITE POSTERIOR causes an allergic reaction or oversensitivity, discontinue its use immediately.
5) Use examination gloves (plastic, vinyl, or latex) at all times when handling ESTELITE POSTERIOR to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If ESTELITE POSTERIOR comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.

- 6) Avoid contact of ESTELITE POSTERIOR with eyes, mucosal membrane, skin and clothing.
- If ESTELITE POSTERIOR comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
- If ESTELITE POSTERIOR comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration.
- If ESTELITE POSTERIOR comes into contact with the skin or clothing, immediately saturate the area with alcohol-soaked cotton swab or gauze.
- Instruct the patient to rinse his mouth immediately after treatment.
7) ESTELITE POSTERIOR should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
8) To avoid the unintentional ingestion of ESTELITE POSTERIOR, do not leave it unsupervised within the reach of patients and children.
9) Clean the placement instruments with alcohol after use.
10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn.

■ INDICATIONS FOR FILLING AND CURING

Be sure to light-cure ESTELITE POSTERIOR extraorally and check the time needed for complete hardening of ESTELITE POSTERIOR with your light curing unit before performing the bonding procedure. The following table illustrates the relationship between curing time and increment depth.

Relationship between curing time and increment depth:

Light Type	Intensity (mW/cm ²)	Curing Time (seconds)	Increment Depth (mm) ⁽¹⁾
Halogen	400	10	1.9
		20	2.3
	600	10	2.0
		20	2.4
LED	800	5	1.8
		10	2.4
	260	10	1.9
		20	2.0
Plasma Arc	900	5	1.8
		10	2.0
	950	3	2.1

1) Increment depth was determined on the basis of test results performed in accordance with section 7.10 "depth of cure" of ISO4049.

- STORAGE**
1) Store ESTELITE POSTERIOR at a temperature between 0 - 25°C (32 - 77°F).
2) AVOID direct exposure to light and heat.
3) DO NOT use ESTELITE POSTERIOR after the indicated date of expiration on the syringe or package.

■ DISPOSAL
To safely dispose of excess ESTELITE POSTERIOR, extrude unused portion from SYRINGE and light-cure before disposal.

■ CLINICAL PROCEDURE

1. Cleaning
Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinses with water.
2. Shade Selection
Select the appropriate shade of ESTELITE POSTERIOR prior to placing the rubber dam.
- Complete the shade selection within 5 minutes; desiccated teeth are lighter than wet or moist teeth.
- Lightness (color value) is primary importance on shade selection.
- When wavering in shade selection, select the lighter (high value) shade.
- In case of whitened teeth select the shade a few weeks after the completion of whitening. Whitened teeth may become slight darker.

3. Isolation

A rubber dam is the preferred method to isolate the tooth from contamination.

4. Cavity Preparation

Prepare the cavity and rinse with water. Add chamfers to the margins of posterior preparations (class I, II) because chamfers assist in diminishing demarcations between the cavity margins and restoration, thereby enhancing both esthetics and retention.

- In case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent according to manufacturer's instructions.

5. Pulp Protection

Glass ionomer lining or calcium hydroxide should be applied if the cavity is in close proximity to the pulp. DO NOT USE EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of bonding agent and/or ESTELITE POSTERIOR.

6. Bonding system

Apply light-cured bonding system according to manufacturer's instructions. We recommend the use of TOKUYAMA BOND FORCE which is a one component, self-etching, light-cured dental adhesive.

- Do not use self-cured bonding systems. We advise that light-cured bonding systems should be used with ESTELITE POSTERIOR since other types of bonding systems such as self-cured or dual-cured bonding systems may not provide effective adhesion when used with light-cured composite resins. If you want to use self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.

7. Dispensing

Remove the syringe cap. Extrude the paste onto a mixing pad by turning the handle clockwise. Re-cap the syringe immediately after dispensing.

- DO NOT apply unnecessary force to the syringe immediately after removing from the refrigerator.

8. Filling and contouring

Fill the cavity incrementally using a suitable placing instrument. Increments should not exceed the indicated depth (please refer to the aforementioned table).

- DO NOT mix ESTELITE POSTERIOR with other brands of resin composite or other shades of ESTELITE POSTERIOR to avoid incomplete cure or entrapment of air bubbles.

9. Curing

Light-cure each increment for at least the indicated time (please see aforementioned table), keeping the curing light tip within a distance of 2 mm from the increment.

- If other brands of composite resins are layered over the cured composite, follow their instructions.

10. Finishing

Shape and polish the restoration. For finishing, use fine finishing diamond points. To diminish the demarcation between the composite and the natural tooth, rotate a 12 fluted carbide bur slowly along the margins without water. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tool, such as PoGo™ (DENTSPLY/Caulk), Sof-Lex™ (3M ESPE), Identoflex® HiLuster Dia Polishers (Kerr), or D-FINE™ Hybrid Diamond (COSMEDENT, Inc.).
* PoGo, Sof-Lex, Identoflex and D-FINE are not registered trademarks of Tokuyama Dental Corporation.

11. Curing

Light-cure each increment for at least the indicated time (please see aforementioned table), keeping the curing light tip within a distance of 2 mm from the increment.

- If other brands of composite resins are layered over the cured composite, follow their instructions.

12. Finishing

Wet shape the composite resin with water and immediately contact an ophthalmologist.

- If ESTELITE POSTERIOR comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration.

- If ESTELITE POSTERIOR comes into contact with the skin or clothing, immediately saturate the area with alcohol-soaked cotton swab or gauze.

- Instruct the patient to rinse his mouth immediately after treatment.

7) ESTELITE POSTERIOR should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.

8) To avoid the unintentional ingestion of ESTELITE POSTERIOR, do not leave it unsupervised within the reach of patients and children.

9) Clean the placement instruments with alcohol after use.

10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn.

■ INDICATIONS FOR FILLING AND CURING

Be sure to light-cure ESTELITE POSTERIOR extraorally and check the time needed for complete hardening of ESTELITE POSTERIOR with your light curing unit before performing the bonding procedure. The following table illustrates the relationship between curing time and increment depth.

Relationship between curing time and increment depth:

Light Type	Intensity (mW/cm ²)	Curing Time (seconds)	Increment Depth (mm) ⁽¹⁾
Halogen	400	10	1.9
		20	2.3
	600	10	2.0
		20	2.4
LED	800	5	1.8
		10	2.4
	260	10	1.9
		20	2.0
Plasma Arc	900	5	1.8
		10	2.0
	950	3	2.1

1) Increment depth was determined on the basis of test results performed in accordance with section 7.10 "depth of cure" of ISO4049.

- STORAGE**
1) Store ESTELITE POSTERIOR at a temperature between 0 - 25°C (32 - 77°F).
2) AVOID direct exposure to light and heat.
3) DO NOT use ESTELITE POSTERIOR after the indicated date of expiration on the syringe or package.

■ DISPOSAL

To safely dispose of excess ESTELITE POSTERIOR, extrude unused portion from SYRINGE and light-cure before disposal.

■ CLINICAL PROCEDURE

1. Cleaning
Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinses with water.
2. Shade Selection
Select the appropriate shade of ESTELITE POSTERIOR prior to placing the rubber dam.
- Complete the shade selection within 5 minutes; desiccated teeth are lighter than wet or moist teeth.
- Lightness (color value) is primary importance on shade selection.
- When wavering in shade selection, select the lighter (high value) shade.
- In case of whitened teeth select the shade a few weeks after the completion of whitening. Whitened teeth may become slight darker.

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