SILICONE REMOVER

Last updated 2019-10-31 Rev.A

ENGLISH

•••••

READ ALL THE INSTRUCTIONS BEFORE USING SILICONE REMOVER.

INSTRUCTIONS FOR USE

SILICONE REMOVER is an auxiliary agent used to help to remove silicone-based resilient denture liner from the acrylic denture and to remove silicone rubber impression material from the impression resin tray. SILICONE REMOVER does not inhibit the curing of silicone-based material in reline material and the reimpression.

PRECAUTIONS

- SILICONE REMOVER should only be used by licensed dental care professionals; dentists, dental hygienists or dental technicians.
- Use SILICONE REMOVER only as instructed herein. Do not use SILICONE REMOVER for any application other than those noted in these instructions.
- Carefully read and follow all manufacturer’s instructions prior to using this product.

■Contents

Silicone Remover (7 mL)

■Clinical Application

- Removal of silicone-based resilient denture liner, such as SOFRELINER TOUGH, from dentures.
- Removal of silicone rubber impression materials from the impression resin trays.

■Safety Precautions

- The use of examination gloves (plastic, vinyl, latex) and protective glasses is recommended at all times to avoid the possibility of allergic reactions from acrylic-based monomers. If SILICONE REMOVER comes in contact with the examination gloves, immediately dispose of the examination gloves and wash hands thoroughly with water. In case of contact with skin, wash the skin with soap, and then thoroughly flush with water.
- Ensure that SILICONE REMOVER does not come in contact with patient. In case of contact with the patient’s skin, wash the skin with soap, and then thoroughly flush with water
- AVOID CONTACT WITH EYES. In case of contact with eyes, thoroughly flush eyes with water, and immediately contact a physician for assistance.
- Do not apply SILICONE REMOVER to the patient’s mouth or gums. SILICONE REMOVER must be applied only to the denture or impression resin tray and always with the denture or impression resin tray outside of the patient’s mouth.
- SILICONE REMOVER is a highly flammable and volatile liquid. Do not use near an open flame. Use it in a well-ventilated area.
- To avoid the risk of cross-infection, do not contact the nozzle of bottle with the denture or relining material.

■Clinical Precautions

- SILICONE REMOVER may spatter during the peeling procedure. In case of contact with eyes, thoroughly flush eyes with water, and immediately contact a physician for assistance. In case of contact with the patient’s skin, wash the skin with soap, and then thoroughly flush with water
- Always use a clean, new plastic cup and brush when using SILICONE REMOVER. Residual dental materials in previously used cups and brushes, such as self curing resin, contain residual catalysts, and may cure SILICONE REMOVER..
- Do not use this product if the safety seal is broken or if the contents have been tampered with.
- Do not immerse the denture or resin tray into SILICONE REMOVER liquid. Otherwise, the surface of the acrylic resin will be unusable.
- Do not use the liquid of self-curing resin instead of SILICONE REMOVER. Additives in the liquid may influence the adhesion to the denture on relining.
- SILICONE REMOVER may not be used for the removal of acrylic resin based soft denture liner and tissue conditioner.

■Storage Instructions


- Store the material at 0-25°C (32-77°F) at all times.
- DO NOT use the material past its expiration date.
- Replace caps tightly when finished.
- To avoid accidental ingestion or spillage, ALWAYS keep out of the reach of children and patients.

■Disposal Instructions

Any remaining SILICONE REMOVER should be absorbed with gauze or cotton, and disposed of in accordance with Federal and/or local regulations. Avoid any inhalation of the volatile substance.

■Directions For Use

- Expose the acrylic surface by cutting the silicon material 2-5 mm in width and 10-20 mm in length with a surgical knife (see Fig. 1).
- Put a few drops of SILICONE REMOVER in the formed concave portion and allow it to sit for 5 – 10 seconds (see Fig. 2).
- Add a few more drops of SILICONE REMOVER, put the instrument or spatula in the boundary between silicone and acrylic resin and then remove the silicone material from the acrylic resin with the instrument or spatula (see Fig. 3). If necessary, add a few more drops of SILICONE REMOVER to the formed concave portion to aid in removal.

Relining and Reimpression:

Relining or reimpression is performed as follows;

- Remove the surface of adhesion area of the denture or resin tray by using a grinder.
- Allow it to rinse and dry.

The manufacture is not responsible for damage or injury caused by improper use of this product. It is the responsibility of the user to ensure the product is suitable for its intended application.

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