## Resin-based Dental Restorative Material

# PALFIQUE ASTERIA

## ENGLISH

Read all information, precautions and notes before using.

### ■ PRODUCT DESCRIPTION AND GENERAL INFORMATION

1) PALFIQUE ASTERIA is a light-cured, radiopaque composite resin for use in anterior and posterior restorations. It is indicated for all carious lesion classes, including minimal or no-preparation type restorations. PALFIQUE ASTERIA contains 82% by weight (71% by volume) silica-zirconia filler and composite filler. A high filler load offers decreased polymerization shrinkage. Every inorganic filler contained in PALFIQUE ASTERIA is a spherical filler (mean particle size: 200 nm, particle size range: 100 to 300 nm) that enables excellent gloss retention and wear resistance. PALFIQUE ASTERIA contains Bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), Bisphenol A polyethoxy methacrylate (Bis-MPEPP), 1,6-bis(methacryl-ethyloxy carbonylamino)trimethyl hexane (UDMA), Triethylene glycol dimethacrylate (TEGDMA), Mequinol, Dibutyl hydroxyl toluene, and UV absorber.

2) PALFIQUE ASTERIA incorporates Radical-Amplified Photopolymerization initiator technology (RAP technology). RAP technology facilitates a shortened light curing time and ample working time (90 sec.). Please see the table depicting the relationship between curing time and increment depth (■ INDICATIONS FOR FILLING AND CURING).

3) PALFIQUE ASTERIA is provided in SYRINGE.

### ■ SHADE

PALFIQUE ASTERIA contains 7 Body shades and 5 Enamel shades. The PALFIQUE ASTERIA shading system is designed to mimic a wide shade range of the natural dentition by layering the appropriate shades selected from 2 categories: Body and Enamel.

Component	Shade	Tip
Body	A1B, A2B, A3B, A3.5B, A4B, B3B, BL	Body shades are designed for reconstructing the dentin layer. The Body shades should cover all enamel margins except the incisal area.
	Natural Enamel (NE)	NE is recommended to restore translucency in the incisal area in most cases.
	White Enamel (WE)	WE is recommended for the proximal wall. WE is suggested as an alternative to NE in whiter cases.
Enamel	Yellow Enamel (YE)	YE is designed to mimic discolored enamel.
	Trans Enamel (TE)	TE is the most translucent in PALFIQUE ASTERIA. This shade is suggested as an alternative to NE in highly translucent cases.
	Occlusal Enamel (OcE)	OcE is recommended for the occlusal surface. OcE has exceptional sculptability to shape occlusal cusps and fissures.

### ■ INDICATIONS

- Direct anterior and posterior restorations including occlusal surfaces

- Direct bonded composite veneer

- Diastema closure

- Repair of porcelain/composite

### ■ CONTRAINDICATIONS

PALFIQUE ASTERIA contains methacrylic monomers. DO NOT use PALFIQUE ASTERIA for patients allergic to or hypersensitive to methacrylic and related monomers or to any of the other ingredients.

### ■ PRECAUTIONS

1) DO NOT use PALFIQUE ASTERIA for any purposes other than those listed in these instructions. Use PALFIQUE ASTERIA only as directed herein.

2) PALFIQUE ASTERIA is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor use by non-dental care professionals.

3) DO NOT use PALFIQUE ASTERIA if the safety seals are broken or appear to have been tampered with.

4) If PALFIQUE ASTERIA causes an allergic reaction or hypersensitivity, discontinue its use immediately.

5) Use examination gloves (plastic, vinyl or latex) at all times when handling PALFIQUE ASTERIA to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If PALFIQUE ASTERIA comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.

6) Avoid contact of PALFIQUE ASTERIA with eyes, mucosal membrane, skin and clothing.

- If PALFIQUE ASTERIA comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.

- If PALFIQUE ASTERIA comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.

- If PALFIQUE ASTERIA comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.

- Instruct the patient to rinse his mouth immediately after treatment.

7) PALFIQUE ASTERIA should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.

8) To avoid the unintentional ingestion of PALFIQUE ASTERIA, do not leave it unsupervised within the reach of patients and children.

9) Clean the placement instruments and brushes with alcohol after use.

10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.

11) This product is designed to be used at room temperature (18 - 30°C / 62 - 84°F). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

### ■ INDICATIONS FOR FILLING AND CURING

PALFIQUE ASTERIA is designed to be cured by either a halogen, LED or Plasma Arc light with a wavelength of 400-500 nm. Be sure to light-cure PALFIQUE ASTERIA extra-orally and check the time needed for complete hardening of PALFIQUE ASTERIA with your light-curing unit before performing the bonding procedure. The following table summarizes the relationship between curing time and increment depth.

Relationship between curing time and increment depth:

Light type	Intensity (mW/cm²)	Curing time (Seconds)	Increment depth (mm) <sup>1)</sup>	
			A1B, A2B, A3B, A3.5B, B3B, BL	A4B
Halogen	400	10	1.7	— <sup>2)</sup>
		20	2.1	1.8
	600	10	2.0	1.7
		20	2.4	2.0
	800	10	2.1	1.8
		20	2.5	2.1
LED	400	10	1.7	— <sup>2)</sup>
		20	2.0	1.7
	600	10	1.8	1.6
		20	2.2	1.9
	1000	10	2.0	1.8
Plasma Arc	950	6	1.9	1.8

1) Increment depth was determined on the basis of test results performed in accordance with section 7.10 "depth of cure" of ISO4049.

2) Curing time must be at least 20 seconds.

### ■ STORAGE

1) Store PALFIQUE ASTERIA at temperatures between 0 - 25°C (32 - 77°F).

2) AVOID direct exposure to light and heat.

3) DO NOT use PALFIQUE ASTERIA after the indicated date of expiration on the syringe package.

### ■ DISPOSAL

To safely dispose of excess PALFIQUE ASTERIA, extrude unused portion from SYRINGE and light-cure before disposal.

### ■ CUSTOM SHADE GUIDE

See the instructions for PALFIQUE ASTERIA Custom Shade Guide prior use.

### ■ CLINICAL PROCEDURE

#### 1. Cleaning

Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste, and then rinse with water.

#### 2. Shade Selection

Select the appropriate shade of PALFIQUE ASTERIA using the custom shade guide prepared in advance.

- Complete the shade selection within 5 minutes; teeth become whiter when dehydrated.

- Lightness (color value) is the primary importance on shade selection.

- In the case of whitened teeth, select the shade a few weeks later after the completion of whitening.

Whitened teeth tend to become slightly darker with the elapse of time.

- Please see the table depicting the shade and tip (■ SHADE).

#### 3. Isolation

A rubber dam is the preferred method of isolation.

#### 4. Cavity Preparation

Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), bevels assist in erasing demarcations between the cavity margins and the restoration, thereby enhancing both esthetics and retention.

- Scalloped (undulant) bevel may be desirable depending on the esthetics and retention.

- In the case where no cavity preparation has been made (caries-free cervical defects), clean the tooth surface with a rubber cup and a fluoride-free cleaning paste. Rinse thoroughly with water.

- In the case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water, air dry thoroughly and treat with a silane coupling agent according to manufacturer's instructions.

#### 5. Pulp Protection

Protect the pulp appropriately, if the cavity is in close proximity to the pulp. In the case of pulp exposure, apply calcium hydroxide for pulp capping. DO NOT USE EUGENOL-BASED MATERIALS to protect the pulp as these materials will inhibit curing of PALFIQUE ASTERIA.

#### 6. Bonding System

Apply a bonding system according to its manufacturer's instructions.

#### 7. Dispensing

Remove the syringe cap. Extrude the paste onto a mixing pad by turning the handle clockwise. After dispensing, turn the handle counterclockwise by a half to full turn to release residual pressure inside the syringe and re-cap the syringe immediately.

- DO NOT apply unnecessary force to the syringe immediately after removing from the refrigerator.

#### 8. Filling and contouring

Fill the cavity incrementally. Increments should not exceed the indicated curing depth (please refer to the aforementioned table).

- In case that characterization is required, tints (such as ESTELITE COLOR) can be used. Apply the tint according to manufacturer's instructions.

- DO NOT mix PALFIQUE ASTERIA with other brands of composite resin to avoid incomplete cure or entrainment of air bubbles. In order to avoid air bubbles entrapment, DO NOT mix PALFIQUE ASTERIA with another shade of the paste.

#### 9. Curing

Light-cure each increment for at least the indicated time (please see aforementioned table), keeping the curing light tip within a distance of 2 mm from the increment.

- If other brands of composite resins are layered over the cured composite, follow its manufacturer's instructions.

#### 10. Finishing

Shape and polish the restoration. For finishing, use finishing discs and/or fine finishing diamond points. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tools. For final polishing, use felt discs or cotton wheels with polishing paste, or suitable polishing tools.

### IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the dental professional to ensure the product is suitable for application before use.

Specifications are subject to change without notice. When the product specification changes, the instructions and precautions may change also.

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