

# TOKUYAMA UNIVERSAL PRIMER

Comprehensive Solution

Single Source  
Solution

Versatile


Cost Effective




# Everything from A-Z

TOKUYAMA UNIVERSAL PRIMER allows you to use it for pretreatment of all materials, such as;

- Bonding of inlays, onlays, crowns, bridges, veneers, cores, artificial teeth, orthodontic brackets etc.
- Repair of ceramic, metal and composite materials (including porcelain fused metals).
- Bonding of denture resin to metal base, clasp or attachment.
- Repair of denture with metal base, clasp or attachment.
- Bonding of opaque resin to a metal base in the fabrication of resin-faced stainless steel crowns.




## Simple Steps



## Reliable


TOKUYAMA UNIVERSAL PRIMER shows high bond strength to a variety of materials



\*\* Source: Tokuyama Dental R&D

\* Not a registered trademark of Tokuyama Dental Corp.

TOKUYAMA UNIVERSAL PRIMER shows high bond strength to all prosthetic materials



## Product Package Available

### TOKUYAMA UNIVERSAL PRIMER KIT

- PRIMER A 2mL
- PRIMER B 2mL



Specifications are subject to change without obligation on the part of the Manufacturer.

Manufactured by  
**Tokuyama Dental Corporation**

38-9, Taitou 1-chome, Taitou-ku, Tokyo, Japan  
TEL: +81-3-3835-2261 FAX: +81-3-3835-2265  
URL: <http://www.tokuyama-dental.com>

