

# SAFETY DATA SHEET

## TOKUYAMA EE-BOND, EE-BOND

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

**Product name** TOKUYAMA EE-BOND, EE-BOND

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

**Identified uses** 【Medical device】 Dental adhesive. For dental professionals only.

**Uses Not Recommended** Use only for intended applications.

#### 1.3. Details of the supplier of the safety data sheet

**Supplier** Tokuyama Dental Italy S.r.l.  
Via Chizzalunga, 1,  
36066 Sandrigo, Vicenza, Italy  
TEL: +39-0444-659650  
FAX: +39-0444-750345

**Contact person** <http://www.tokuyama-dental.com/tdc/contact.html>

**Manufacturer** Tokuyama Dental Corporation  
38-9, Taitou 1-chome, Taitou-ku, Tokyo  
110-0016, Japan  
TEL: +81-3-3835-2261  
FAX: +81-3-3835-2265

#### 1.4. Emergency telephone number

**Emergency telephone** Regional Medicines and Poisons Information Centre NI,  
Pharmacy Department  
Royal Hospital Suite, Grosvenor Road, Belfast  
Emergency telephone: 844 892 0111  
Telephone: +44 28 90 63 2032  
Fax: +44 28 90 24 8030  
E-mail address: [nirdic.nirdic@belfasttrust.hscni.net](mailto:nirdic.nirdic@belfasttrust.hscni.net)

### SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture

##### Classification (EC 1272/2008)


**Physical hazards** Flam. Liq. 2 - H225

**Health hazards** Skin Irrit. 2 - H315: Eye Irrit. 2 - H319: Skin Sens. 1 - H317: STOT SE 3 - H336

**Environmental hazards** Not Classified

#### 2.2. Label elements

##### Pictogram



**Signal word** Danger

**Hazard statements** H225 Highly flammable liquid and vapour.  
H315 Causes skin irritation.  
H317 May cause an allergic skin reaction.  
H319 Causes serious eye irritation.  
H336 May cause drowsiness or dizziness.

## TOKUYAMA EE-BOND, EE-BOND

**Precautionary statements** P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.  
 P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.  
 P403+P235 Store in a well-ventilated place. Keep cool.  
 P501 Dispose of contents/ container in accordance with local regulations.

**Contains** (1-METHYLETHYLIDENE)BIS[4,1-PHENYLENEOXY(2-HYDROXY-3,1-PROPANEDIYL)]  
 BISMETHACRYLATE, 2-HYDROXYETHYL METHACRYLATE,  
 2,2'-ETHYLENEDIOXYDIETHYL DIMETHACRYLATE

**Supplementary precautionary statements** P240 Ground/ bond container and receiving equipment.  
 P241 Use explosion-proof electrical equipment.  
 P242 Use only non-sparking tools.  
 P243 Take precautionary measures against static discharge.  
 P261 Avoid breathing vapour/ spray.  
 P264 Wash contaminated skin thoroughly after handling.  
 P271 Use only outdoors or in a well-ventilated area.  
 P272 Contaminated work clothing should not be allowed out of the workplace.  
 P302+P352 IF ON SKIN: Wash with plenty of water.  
 P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/ shower.  
 P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.  
 P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
 P312 Call a POISON CENTER/ doctor if you feel unwell.  
 P321 Specific treatment (see medical advice on this label).  
 P332+P313 If skin irritation occurs: Get medical advice/ attention.  
 P333+P313 If skin irritation or rash occurs: Get medical advice/ attention.  
 P337+P313 If eye irritation persists: Get medical advice/ attention.  
 P362+P364 Take off contaminated clothing and wash it before reuse.  
 P370+P378 In case of fire: Use foam, carbon dioxide, dry powder or water fog to extinguish.  
 P403+P233 Store in a well-ventilated place. Keep container tightly closed.

**2.3. Other hazards** This product does not contain any substances classified as PBT or vPvB.  
 Endocrine disrupting properties: Not known.

<b>SECTION 3: Composition/information on ingredients</b>
--

**3.2. Mixtures**

<b>PROPAN-2-OL</b>	<b>30-60%</b>
CAS number: 67-63-0	EC number: 200-661-7
<b>Classification</b>	
Flam. Liq. 2 - H225	
Eye Irrit. 2 - H319	
STOT SE 3 - H336	

## TOKUYAMA EE-BOND, EE-BOND

<b>(1-METHYLETHYLIDENE)BIS[4,1-PHENYLENEOXY(2-HYDROXY-3,1-PROPANEDIYL)] BISMETHACRYLATE</b> <span style="float: right;"><b>10-30%</b></span>
CAS number: 1565-94-2                      EC number: 216-367-7
<b>Classification</b> Skin Sens. 1 - H317 Aquatic Chronic 3 - H412
<b>2-PROPENOIC ACID, 2-METHYL-, 2-HYDROXYETHYL ESTER, PHOSPHATE</b> <span style="float: right;"><b>10-30%</b></span>
CAS number: 52628-03-2                      EC number: 258-053-2
<b>Classification</b> Skin Irrit. 2 - H315 Eye Irrit. 2 - H319
<b>2-HYDROXYETHYL METHACRYLATE</b> <span style="float: right;"><b>10-30%</b></span>
CAS number: 868-77-9                      EC number: 212-782-2
<b>Classification</b> Skin Irrit. 2 - H315 Eye Irrit. 2 - H319 Skin Sens. 1 - H317
<b>2,2'-ETHYLENEDIOXYDIETHYL DIMETHACRYLATE</b> <span style="float: right;"><b>5-10%</b></span>
CAS number: 109-16-0                      EC number: 203-652-6
<b>Classification</b> Skin Irrit. 2 - H315 Eye Irrit. 2 - H319 Skin Sens. 1 - H317
<b>DL-BORNANE-2,3-DIONE</b> <span style="float: right;"><b>&lt;1%</b></span>
CAS number: 10373-78-1                      EC number: 233-814-1
<b>Classification</b> Skin Irrit. 2 - H315 Eye Irrit. 2 - H319 STOT SE 3 - H335
<b>2,6-DI-tert-BUTYL-p-CRESOL</b> <span style="float: right;"><b>&lt;1%</b></span>
CAS number: 128-37-0                      EC number: 204-881-4 M factor (Chronic) = 1
<b>Classification</b> Aquatic Chronic 1 - H410

# TOKUYAMA EE-BOND, EE-BOND

<b>DIPHENYL(2,4,6-TRIMETHYLBENZOYL)PHOSPHINE OXIDE</b>	<b>&lt;1%</b>
CAS number: 75980-60-8	EC number: 278-355-8

<b>Classification</b> Repr. 2 - H361f
--

<b>MEQUINOL</b>	<b>&lt;1%</b>
CAS number: 150-76-5	EC number: 205-769-8

<b>Classification</b> Acute Tox. 4 - H302 Eye Irrit. 2 - H319 Skin Sens. 1 - H317
--

The full text for all hazard statements is displayed in Section 16.

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

<b>General information</b>	Move affected person to fresh air at once. Get medical attention if any discomfort continues.
<b>Inhalation</b>	Move affected person to fresh air at once.
<b>Ingestion</b>	Rinse mouth thoroughly with water. Give plenty of water to drink. Get medical attention if any discomfort continues.
<b>Skin contact</b>	Wash skin thoroughly with soap and water. Take off contaminated clothing and wash it before reuse.
<b>Eye contact</b>	Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Get medical attention if irritation persists after washing.

### 4.2. Most important symptoms and effects, both acute and delayed

<b>Inhalation</b>	Vapours may cause drowsiness and dizziness.
<b>Skin contact</b>	Skin irritation. Allergic rash.
<b>Eye contact</b>	May cause severe eye irritation.

### 4.3. Indication of any immediate medical attention and special treatment needed

<b>Notes for the doctor</b>	No specific recommendations.
-----------------------------	------------------------------

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

<b>Suitable extinguishing media</b>	Extinguish with the following media: Foam. Carbon dioxide or dry powder.
-------------------------------------	--

### 5.2. Special hazards arising from the substance or mixture

<b>Specific hazards</b>	Irritating gases or vapours. The product is highly flammable.
<b>Hazardous combustion products</b>	No known hazardous decomposition products.

### 5.3. Advice for firefighters

<b>Protective actions during firefighting</b>	Avoid breathing fire gases or vapours.
---	--

# TOKUYAMA EE-BOND, EE-BOND

**Special protective equipment for firefighters** Wear chemical protective suit.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

**Personal precautions** Wear protective clothing as described in Section 8 of this safety data sheet. No smoking, sparks, flames or other sources of ignition near spillage.

### 6.2. Environmental precautions

**Environmental precautions** Do not discharge into drains or watercourses or onto the ground.

### 6.3. Methods and material for containment and cleaning up

**Methods for cleaning up** Eliminate all sources of ignition. No smoking, sparks, flames or other sources of ignition near spillage. Provide adequate ventilation. Absorb in vermiculite, dry sand or earth and place into containers. Use only non-sparking tools. Ground/bond container and receiving equipment. Avoid the spillage or runoff entering drains, sewers or watercourses.

### 6.4. Reference to other sections

**Reference to other sections** For waste disposal, see Section 13.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

**Usage precautions** Wear protective clothing as described in Section 8 of this safety data sheet. All handling should only take place in well-ventilated areas. Avoid inhalation of vapours/spray and contact with skin and eyes. Take precautionary measures against static discharges. Use explosion-proof electrical, ventilating and lighting equipment.

### 7.2. Conditions for safe storage, including any incompatibilities

**Storage precautions** Store in a cool and well-ventilated place. Protect from sunlight. Store at temperatures between 0°C/10°F and 32°C/50°F.

**Storage class** Flammable liquid storage.

### 7.3. Specific end use(s)

**Specific end use(s)** The identified uses for this product are detailed in Section 1.2.

## SECTION 8: Exposure Controls/personal protection

### 8.1. Control parameters

#### Occupational exposure limits

##### PROPAN-2-OL

Long-term exposure limit (8-hour TWA): WEL 400 ppm 999 mg/m<sup>3</sup>

Short-term exposure limit (15-minute): WEL 500 ppm 1250 mg/m<sup>3</sup>

##### 2,6-DI-tert-BUTYL-p-CRESOL

Long-term exposure limit (8-hour TWA): WEL 10 mg/m<sup>3</sup>

WEL = Workplace Exposure Limit

### 8.2. Exposure controls

**Appropriate engineering controls** Provide adequate general and local exhaust ventilation. Use explosion-proof electrical, ventilating and lighting equipment.

**Eye/face protection** The following protection should be worn: Chemical splash goggles.

**Hand protection** Wear protective gloves.

## TOKUYAMA EE-BOND, EE-BOND

<b>Other skin and body protection</b>	Wear appropriate clothing to prevent any possibility of skin contact.
<b>Hygiene measures</b>	Wash hands at the end of each work shift and before eating, smoking and using the toilet. Wash promptly with soap and water if skin becomes contaminated. Promptly remove any clothing that becomes contaminated. Wash contaminated clothing before reuse. Use appropriate skin cream to prevent drying of skin. When using do not eat, drink or smoke.
<b>Respiratory protection</b>	No specific recommendations. Respiratory protection must be used if the airborne contamination exceeds the recommended occupational exposure limit.

### SECTION 9: Physical and Chemical Properties

#### 9.1. Information on basic physical and chemical properties

<b>Appearance</b>	Liquid
<b>Colour</b>	Light (or pale) yellow
<b>Odour</b>	Characteristic
<b>Odour threshold</b>	Not available.
<b>pH</b>	pH (diluted solution): 2.3
<b>Melting point</b>	Not available.
<b>Initial boiling point and range</b>	Not available.
<b>Flash point</b>	20.0°C
<b>Evaporation rate</b>	Not available.
<b>Evaporation factor</b>	Not available.
<b>Flammability (solid, gas)</b>	Not available.
<b>Upper/lower flammability or explosive limits</b>	Not applicable.
<b>Other flammability</b>	Not available.
<b>Vapour pressure</b>	Not available.
<b>Vapour density</b>	Not available.
<b>Relative density</b>	1.00
<b>Solubility(ies)</b>	Not available.
<b>Partition coefficient</b>	Not available.
<b>Auto-ignition temperature</b>	Not available.
<b>Decomposition Temperature</b>	Not available.
<b>Viscosity</b>	Not available.
<b>Explosive properties</b>	Not explosive.
<b>Oxidising properties</b>	Not available.

#### 9.2. Other information

<b>Other information</b>	Not available.
--------------------------	----------------

### SECTION 10: Stability and reactivity

#### 10.1. Reactivity

## TOKUYAMA EE-BOND, EE-BOND

**Reactivity** There are no known reactivity hazards associated with this product. See the other subsections of this section for further details.

### 10.2. Chemical stability

**Stability** The product is highly flammable.

### 10.3. Possibility of hazardous reactions

**Possibility of hazardous reactions** Not known.

### 10.4. Conditions to avoid

**Conditions to avoid** Keep away from heat, sparks and open flame. Protect against direct sunlight.

### 10.5. Incompatible materials

**Materials to avoid** Strong oxidising agents.

### 10.6. Hazardous decomposition products

**Hazardous decomposition products** Heating may generate the following products: Toxic gases or vapours.

## SECTION 11: Toxicological information

### 11.1 Information on toxicological effects

#### **Acute toxicity - Ingestion**

Calculation method : Not classified.

#### **Acute toxicity - Skin Contact**

Calculation method : Not classified. No data.

#### **Acute toxicity - Inhalation**

Calculation method : Not classified. No data.

#### **Skin corrosion/irritation**

Calculation method : Causes skin irritation.

#### **Serious eye damage/irritation**

Calculation method : Causes serious eye irritation.

#### **Skin sensitization**

Calculation method : May cause an allergic skin reaction.

#### **Respiratory sensitization**

Calculation method : Not classified.

#### **Germ cell mutagenicity**

Calculation method : Not classified.

#### **Carcinogenicity**

Calculation method : Not classified.

#### **Reproductive toxicity**

Calculation method : Not classified.

#### **Lactation**

Calculation method : Not classified.

#### **STOT - single exposure**

Calculation method : May cause drowsiness or dizziness.

# TOKUYAMA EE-BOND, EE-BOND

## **STOT - repeated exposure**

Calculation method : Not classified.

## **Aspiration hazard**

Calculation method : Not classified. Not applicable.

## **11.2 Other information**

Route(s) of entry Skin, Eyes, Ingestion, Inhalation, Health effects: See section 4.2.

Endocrine disrupting properties: Not known.

## **SECTION 12: Ecological Information**

### **12.1 Toxicity**

**Toxicity - Aquatic invertebrates** Not classified.

**Toxicity - Fish** Not classified.

**Toxicity - Algae** Not classified.

**Toxicity - Sediment Compartment** Not classified.

**Toxicity - Terrestrial Compartment** Not classified.

### **12.2 Persistence and degradability**

Not known.

### **12.3 Bioaccumulative potential**

Not known.

### **12.4 Mobility in soil**

Not known.

### **12.5 Results of PBT and vPvB assessment**

This product does not contain any substances classified as PBT or vPvB.

### **12.6 Endocrine disrupting properties**

Not known.

### **12.7 Other adverse effects**

Not known.

## **SECTION 13: Disposal considerations**

### **13.1 Waste treatment methods**

Dispose of contents/ container in accordance with the local/ regional/ national/ international regulations.

Must not be disposed of together with household garbage. Dispose at suitable refuse site. Do not allow to enter drains, sewers or watercourses. Empty containers may retain some product residues and hence be potentially hazardous.

### **13.2 Additional Information**

Disposal should be in accordance with local, state or national legislation.

## **SECTION 14: Transport information**

### **14.1. UN number**

**UN No. (ADR/RID)** 1133

**UN No. (IMDG)** 1133

**UN No. (ICAO)** 1133

### **14.2. UN proper shipping name**

**Proper shipping name (ADR/RID)** ADHESIVES



## TOKUYAMA EE-BOND, EE-BOND

Proper shipping name (IMDG) ADHESIVES

Proper shipping name (ICAO) ADHESIVES

Proper shipping name (ADN) ADHESIVES

### 14.3. Transport hazard class(es)

ADR/RID class 3

ADR/RID label 3

IMDG class 3

ICAO class/division 3

Transport labels



### 14.4. Packing group

ADR/RID packing group II

IMDG packing group II

ICAO packing group II

### 14.5. Environmental hazards

### 14.6. Special precautions for user

EmS F-E, S-D

Emergency Action Code •3YE

Hazard Identification Number 33  
(ADR/RID)

Tunnel restriction code (D/E)

### 14.7. Transport in bulk according to Annex II of MARPOL and the IBC Code

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code Not applicable.

Annex II of MARPOL 73/78 and the IBC Code

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

EU legislation (EC) No 1907/2006 (REACH).  
(EC) No 1272/2008 (CLP).  
(EU) No 2015/830.

### 15.2. Chemical safety assessment

No chemical safety assessment has been carried out.

## SECTION 16: Other information

Classification procedures according to Regulation (EC) 1272/2008 Flam. Liq. 2 - H225: On basis of test data. Skin Irrit. 2 - H315: Eye Irrit. 2 - H319: Skin Sens. 1 - H317: STOT SE 3 - H336 : Calculation method.

## TOKUYAMA EE-BOND, EE-BOND

<b>Training advice</b>	Ensure operatives are trained to minimise exposure. Only trained personnel should use this material.
<b>Revision comments</b>	
<b>Revision date</b>	15/12/2022
<b>Revision</b>	4
<b>Hazard statements in full</b>	H225 Highly flammable liquid and vapour. H302 Harmful if swallowed. H315 Causes skin irritation. H317 May cause an allergic skin reaction. H319 Causes serious eye irritation. H335 May cause respiratory irritation. H336 May cause drowsiness or dizziness. H361f Suspected of damaging fertility. H410 Very toxic to aquatic life with long lasting effects. H412 Harmful to aquatic life with long lasting effects.

The information which is contained in this document is based on available data. However, as such has been obtained from various sources, including independent laboratories, it is given without warranty or representation that is complete, accurate and can be relied upon. Tokuyama Dental Corp. has not attempted to conceal in any way the deleterious aspects of the product listed herein, but makes no warranty as to such information.