# SAFETY DATA SHEET ESTELITE COLOR

## 1. Identification

## **Product identifier**

Product name ESTELITE COLOR

## Recommended use of the chemical and restrictions on use

**Restriction on use** Light-cured Characterizing Material. For dental professionals only.

#### Details of the supplier of the safety data sheet

**Supplier** Tokuyama Dental America, Inc.

740 Garden View CT., Suite 200 Encinitas, CA 92024 U.S.A.

Tel: +1-760-942-7211 Fax:+1-760-942-7212

Contact person http://www.tokuyama-dental.com/tdc/contact.htm

Manufacturer Tokuyama Dental Corporation

38-9, Taitou 1-chome, Taitou-ku, Tokyo

110-0016, Japan TEL: +81-3-3835-2261 FAX: +81-3-3835-2265

#### **Emergency telephone number**

Emergency telephone Tokuyama Dental America, Inc. Tel: +1-760-942-7211

Days and hours of operation: 9:00-17:00(PST/PDT), Weekdays, except national holidays in

US.

Language: English only.

#### 2. Hazard identification

#### Classification of the substance or mixture

WHMIS Regulatory Status This product has been classified according to the hazard criteria of the Hazardous Product

Regulations and the SDS contains all required information.

Physical hazards Not Classified

Health hazards Skin Irrit. 2 - H315 Eye Irrit. 2A - H319 Skin Sens. 1 - H317 STOT SE 3 - H335

Environmental hazards Not Classified

## Label elements

#### **Pictogram**



Signal word Warning

Hazard statements H315 Causes skin irritation.

H317 May cause an allergic skin reaction. H319 Causes serious eye irritation. H335 May cause respiratory irritation.

## **ESTELITE COLOR**

**Precautionary statements** P261 Avoid breathing dust.

P264 Wash contaminated skin thoroughly after handling.

P271 Use only outdoors or in a well-ventilated area.

P272 Contaminated work clothing should not be allowed out of the workplace.
P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

P302+P352 IF ON SKIN: Wash with plenty of water.

P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove

contact lenses, if present and easy to do. Continue rinsing.
P312 Call a POISON CENTRE/doctor if you feel unwell.
P321 Specific treatment (see medical advice on this label).
P332+P313 If skin irritation occurs: Get medical advice/ attention.

P333+P313 If skin irritation or rash occurs: Get medical advice/ attention. P337+P313 If eye irritation persists: Get medical advice/ attention.

P362+P364 Take off contaminated clothing and wash it before reuse. P403+P233 Store in a well-ventilated place. Keep container tightly closed.

P405 Store locked up.

P501 Dispose of contents/ container in accordance with local regulations.

## 3. Composition/information on ingredients

#### **Mixtures**

POLY(OXY-1,2-ETHANEDIYL),  $\alpha,\alpha'$ -[(1-METHYLETHYLIDENE)DI-4,1-PHENYLENE]BIS[ $\omega$ -[(2-

METHYL-1-OXO-2-PROPENYL)OXY]-

CAS number: 41637-38-1

2-PROPENOIC ACID, 2-METHYL-, 1,2-ETHANEDIYLBIS(OXY-2,1-ETHANEDIYL) ESTER

CAS number: 109-16-0

2-PROPENOIC ACID, 2-METHYL-, (1-METHYLETHYLIDENE)BIS[4,1-PHENYLENEOXY(2-HYDROXY-3,1-PROPANEDIYL)] ESTER

CAS number: 1565-94-2

TITANIUM DIOXIDE (TiO2)
CAS number: 13463-67-7

2,6-DI-t-BUTYL-4-METHYLPHENOL

CAS number: 128-37-0
PHENOL, 4-METHOXY-

CAS number: 150-76-5

Composition comments

All percentages displayed expressed as weight/weight.

## 4. First-aid measures

## Description of first aid measures

10-30%

5-10%

<1%

<0.1%

<0.1%

< 0.1%

## **ESTELITE COLOR**

General information Move affected person to fresh air at once. Get medical attention if any discomfort continues.

Inhalation Move affected person to fresh air at once. Get medical advice/attention if you feel unwell.

Ingestion Rinse mouth thoroughly with water. Give plenty of water to drink. Try to induce vomiting. Get

medical attention if any discomfort continues.

Skin contact Remove contaminated clothing. Wash skin thoroughly with soap and water. If skin irritation or

rash occurs: Get medical advice/ attention. Wash contaminated clothing before reuse.

Eye contact Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy

to do. Continue rinsing. Get medical attention if irritation persists after washing.

## Most important symptoms and effects, both acute and delayed

Inhalation Irritation of nose, throat and airway.

Ingestion Ingestion may cause severe irritation of the mouth, the esophagus and the gastrointestinal

Skin contact Skin irritation. Allergic rash.

Eye contact May cause severe eye irritation.

#### Indication of any immediate medical attention and special treatment needed

Notes for the doctor Treat symptomatically.

## 5. Fire-fighting measures

#### Extinguishing media

Suitable extinguishing media Extinguish with foam, carbon dioxide, dry powder or water fog.

## Specific hazards arising from the hazardous product

Specific hazards In case of fire, toxic gases may be formed.

Hazardous combustion Thermal decomposition or combustion may liberate carbon oxides and other toxic gases or

products vapours.

Advice for firefighters

Protective actions during

firefighting

Avoid breathing fire gases or vapours.

**Special protective equipment** Wear chemical protective suit.

for firefighters

#### 6. Accidental release measures

#### Personal precautions, protective equipment and emergency procedures

Personal precautions Wear protective clothing as described in Section 8 of this safety data sheet.

**Environmental precautions** 

**Environmental precautions** Avoid the spillage or runoff entering drains, sewers or watercourses.

# Methods and material for containment and cleaning up

Methods for cleaning up Provide adequate ventilation. Collect and place in suitable waste disposal containers and seal

securely.

Reference to other sections

Reference to other sections For waste disposal, see Section 13.

## 7. Handling and storage

## **ESTELITE COLOR**

#### Precautions for safe handling

Usage precautions Wear protective clothing as described in Section 8 of this safety data sheet. Keep away from

heat, sparks and open flame. All handling should only take place in well-ventilated areas.

Avoid inhalation of vapours/spray and contact with skin and eyes.

Advice on general occupational hygiene

Do not eat, drink or smoke when using this product. Wash at the end of each work shift and before eating, smoking and using the toilet. Wash promptly with soap and water if skin becomes contaminated. Promptly remove any clothing that becomes wet or contaminated.

Wash contaminated clothing before reuse.

#### Conditions for safe storage, including any incompatibilities

Storage precautions Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No

smoking. Avoid exposure to high temperatures or direct sunlight. Store in tightly-closed, original container in a dry, cool and well-ventilated place. Store at temperatures between

0°C/32°F and 10°C/50°F.

Specific end use(s)

Specific end use(s) The identified uses for this product are detailed in Section 1.

#### 8. Exposure Controls/personal protection

## **Control parameters**

### Occupational exposure limits

#### **TITANIUM DIOXIDE (TiO2)**

Long-term exposure limit (8-hour TWA): ACGIH 10 mg/m<sup>3</sup>

Α

#### 2,6-DI-t-BUTYL-4-METHYLPHENOL

Long-term exposure limit (8-hour TWA): ACGIH 2 mg/m³ inhalable fraction and vapor

**A**4

#### PHENOL, 4-METHOXY-

Long-term exposure limit (8-hour TWA): ACGIH 5 mg/m³

ACGIH = American Conference of Governmental Industrial Hygienists.

A4 = Not Classifiable as a Human Carcinogen.

#### **Exposure controls**

Appropriate engineering

controls

Provide adequate general and local exhaust ventilation.

**Eve/face protection** The following protection should be worn: Chemical splash goggles.

**Hand protection** Wear protective gloves.

Other skin and body

protection

Wear appropriate clothing to prevent any possibility of skin contact.

Hygiene measures When using do not eat, drink or smoke. Wash at the end of each work shift and before eating,

smoking and using the toilet. Wash promptly with soap and water if skin becomes contaminated. Promptly remove any clothing that becomes contaminated. Wash

contaminated clothing before reuse.

**Respiratory protection**No specific recommendations. Respiratory protection must be used if the airborne

contamination exceeds the recommended occupational exposure limit.

# 9. Physical and Chemical Properties

## Information on basic physical and chemical properties

Appearance Paste

# **ESTELITE COLOR**

Colour Blue, Gray, Lavender, Ochre, Yellow, Red, Pink, White, Yellow white, Dark brown, Clear

Odour Characteristic
Odour threshold Not available.

**pH** Technical impossibility to obtain the data.

Melting pointNot available.Initial boiling point and rangeNot available.Flash pointNot applicable.

Evaporation rate Not available.

Evaporation factor Not available.

Flammability (solid, gas) Not available.

Upper/lower flammability or explosive limits

Not applicable.

Other flammability

Vapour pressure

Not available.

Vapour density

Not available.

Relative density 1.7

Solubility(ies)

Partition coefficient

Not available.

Auto-ignition temperature

Not available.

Decomposition Temperature

Not available.

Viscosity

Not available.

Explosive properties

Not available.

Oxidising properties

Not available.

Other information Not available.

## 10. Stability and reactivity

**Reactivity** Not known. See the other subsections of this section for further details.

**Stability** Stable at normal ambient temperatures and when used as recommended.

Possibility of hazardous

reactions

Not known.

Conditions to avoid Keep away from heat, sparks and open flame. Protect against direct sunlight.

Materials to avoid No specific material or group of materials is likely to react with the product to produce a

hazardous situation.

Hazardous decomposition

products

Heating may generate the following products: Toxic gases or vapours.

# 11. Toxicological information

# **ESTELITE COLOR**

## Information on toxicological effects

Routes of exposure Skin, Eyes, Ingestion, Health effects: See section 4.

Inhalation: Titanium dioxide (TiO2) is kneaded in the product, it does not exist as respirable particles. Due to the physical nature of this product, exposure by this route is unlikely.

Acute toxicity - oral

Notes (oral LD<sub>50</sub>) Not available.

Acute toxicity - dermal

Notes (dermal LD<sub>50</sub>) Not available.

Acute toxicity - inhalation

Notes (inhalation LC<sub>50</sub>) Not available.

Skin corrosion/irritation

**Skin corrosion/irritation** Causes skin irritation.

Animal data Not available.

Serious eye damage/irritation

**Serious eye damage/irritation** Causes serious eye irritation.

Respiratory sensitization

**Respiratory sensitization** Not available.

Skin sensitization

**Skin sensitization** May cause an allergic skin reaction.

Germ cell mutagenicity

Genotoxicity - in vitro Not available.

**Genotoxicity - in vivo** Not available.

Carcinogenicity

Carcinogenicity Based on available data the classification criteria are not met.

**IARC carcinogenicity** Some of the ingredients are listed or exempt.

NTP carcinogenicity Some of the ingredients are listed or exempt.

Reproductive toxicity

Reproductive toxicity - fertility Not available.

Reproductive toxicity -

Not available.

development

Specific target organ toxicity - single exposure

**STOT - single exposure** May cause respiratory irritation.

Specific target organ toxicity - repeated exposure

**STOT - repeated exposure** Not available.

Aspiration hazard

Aspiration hazard Not available.

Toxicological information on ingredients

**TITANIUM DIOXIDE (TiO2)** 

Carcinogenicity

# **ESTELITE COLOR**

IARC carcinogenicity IARC Group 2B Possibly carcinogenic to humans.

**NTP carcinogenicity** Reasonably anticipated to be a human carcinogen.

2,6-DI-t-BUTYL-4-METHYLPHENOL

Carcinogenicity

IARC carcinogenicity IARC Group 3 Not classifiable as to its carcinogenicity to humans.

NTP carcinogenicity Reasonably anticipated to be a human carcinogen.

12. Ecological Information

**Toxicity** No data available.

**Ecological information on ingredients** 

2,6-DI-t-BUTYL-4-METHYLPHENOL

Chronic aquatic toxicity

M factor (chronic) 1

Persistence and degradability

Persistence and degradability No data available.

Bioaccumulative potential

Bioaccumulative potential No data available on bioaccumulation.

Partition coefficient Not available.

Mobility in soil

Mobility No information available.

Other adverse effects

Other adverse effects Not known.

13. Disposal considerations

Waste treatment methods

Disposal methods Dispose of waste to licensed waste disposal site in accordance with the requirements of the

local Waste Disposal Authority.

14. Transport information

General The product is not covered by international regulations on the transport of dangerous goods

(IMDG, IATA, TDG).

UN number

UN No. (International) Not applicable.

UN proper shipping name

Proper shipping name

(International)

Not applicable.

Transport hazard class(es)

Transport Labels

(International)

No transport warning sign required.

## **ESTELITE COLOR**

Packing group

Packing group (International) Not applicable.

**Environmental hazards** 

Environmentally hazardous substance/marine pollutant

No.

Special precautions for user

Not applicable.

Transport in bulk according to Not applicable.

Annex II of MARPOL 73/78

and the IBC Code

#### 15. Regulatory information

Regulatory References According to WHMIS 2015, in compliance with the Hazardous Product Act (HPA, as

amended) and the requirements of the Hazardous Product Regulations (HPR) Hazardous Products Regulations - SOR/2015-17. Transportation of Dangerous Goods Regulations -

SOR/2015-100.

#### 16. Other information

**Training advice** Ensure operatives are trained to minimize exposure. Only trained personnel should use this

material.

**Revision comments** This is the first issue.

Revision date 2018-06-20

Revision 1

Hazard statements in full H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation. H335 May cause respiratory irritation.

The information which is contained in this document is based on available data. However, as such has been obtained from various sources, including independent laboratories, it is given without warranty or representation that is complete, accurate and can be replied upon. Tokuyama Dental Corp. has not attempted to conceal in any way the deleterious aspects of the product listed herein, but makes no warranty as to such information.